



March 1, 2018

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amylor.com

ACCOUNT INFORMATION

P.O. # or Quote: Billing Name:
Contact Name: Shipping Name:
Email Address: Ship Address:
Phone Number: City:
Fax Number: Province:
User Name or Number: Postal Code:

CG Air THERAPEUTIC SEAT CUSHION

14" to 20" WIDE

SINGLE (1) VALVE

- CGA1VSCxxx-x CG Air Single Valve (Includes a Standard 2-way Stretch Covers) \$600
CGA1V3Dxxx-x CG Air Single Valve (Includes a 3D Mesh Covers) \$690
CGA1VICxxx-x CG Air Single Valve (Includes a Incontinence Covers) \$690

DUAL (2) VALVE

- CGA2VSCxxx-x CG Air Dual Valve (Includes a Standard 2-way Stretch Covers) \$600
CGA2V3Dxxx-x CG Air Dual Valve (Includes a 3D Mesh Covers) \$690
CGA2VICxxx-x CG Air Dual Valve (Includes a Incontinence Covers) \$690

SELECT CELL HEIGHT

2" 3" 4"
14" x 14" Qty: 16" x 16" Qty: 18" x 18" Qty: 20" x 18" Qty:
14" x 16" Qty: 16" x 18" Qty: 18" x 20" Qty: 20" x 20" Qty:
16" x 14" Qty: 18" x 16" Qty: 20" x 16" Qty: Other, specify Qty:

22" to 24" WIDE

SINGLE (1) VALVE

- CGA1VSCxxx-x CG Air Single Valve (Includes a Standard 2-way Stretch Covers) \$890
CGA1V3Dxxx-x CG Air Single Valve (Includes a 3D Mesh Covers) \$910
CGA1VICxxx-x CG Air Single Valve (Includes a Incontinence Covers) \$910

DUAL (2) VALVE

- CGA2VSCxxx-x CG Air Dual Valve (Includes a Standard 2-way Stretch Covers) \$890
CGA2V3Dxxx-x CG Air Dual Valve (Includes a 3D Mesh Covers) \$910
CGA2VICxxx-x CG Air Dual Valve (Includes a Incontinence Covers) \$910

SELECT CELL HEIGHT

2" 3" 4"
22" x 18" Qty: 22" x 22" Qty: 24" x 22" Qty: Other, specify Qty:
22" x 20" Qty: 24" x 20" Qty: 24" x 24" Qty:

26" + WIDE

SINGLE (1) VALVE

- CGA1VSCxxx-x CG Air Single Valve (Includes a Standard 2-way Stretch Covers) \$990
CGA1V3Dxxx-x CG Air Single Valve (Includes a 3D Mesh Covers) \$1,010
CGA1VICxxx-x CG Air Single Valve (Includes a Incontinence Covers) \$1,010

DUAL (2) VALVE

- CGA2VSCxxx-x CG Air Dual Valve (Includes a Standard 2-way Stretch Covers) \$990
CGA2V3Dxxx-x CG Air Dual Valve (Includes a 3D Mesh Covers) \$1,010
CGA2VICxxx-x CG Air Dual Valve (Includes a Incontinence Covers) \$1,010

SELECT CELL HEIGHT

2" 3" 4"
26" x 20" Qty: 28" x 20" Qty: 30" x 20" Qty: Other, specify Qty:
26" x 22" Qty: 28" x 22" Qty: 30" x 22" Qty:
26" x 24" Qty: 28" x 24" Qty: 30" x 24" Qty:

REPLACEMENT COVERS

- CGASCxxx-x Standard 2-way Stretch Cover. Specify Dimensions: W " D " H " \$96
CGA3Dxxx-x 3D Mesh Cover. Specify Dimensions: W " D " H " \$144
CGAICxxx-x Incontinence Cover. Specify Dimensions: W " D " H " \$144



CG Air is an Amylor Inc Brand



March 1, 2018

DELUXE

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amylor.com

ACCOUNT INFORMATION

P.O. # or Quote:	_____	Billing Name:	_____
Contact Name:	_____	Shipping Name:	_____
Email Address:	_____	Ship Address:	_____
Phone Number:	_____	City:	_____
Fax Number:	_____	Province:	_____
User Name or Number:	_____	Postal Code:	_____

CG Air DELUXE - THERAPEUTIC CUSTOM SEAT CUSHION

14" to 20" WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX10-20-SC	CG Air Deluxe (Includes a (1) Standard 2-way Stretch Covers) - Below 22" Wide	\$790
<input type="checkbox"/>	CGADLX10-20-3D	CG Air Deluxe (Includes a (1) 3D Mesh Covers) - Below 22" Wide	\$870
<input type="checkbox"/>	CGADLX10-20-IC	CG Air Deluxe (Includes an (1) Incontinence Covers) - Below 22" Wide	\$870

22" to 24" WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX22-24-SC	CG Air Deluxe (Includes a (1) Standard 2-way Stretch Covers) - 22" to 24" Wide	\$890
<input type="checkbox"/>	CGADLX22-24-3D	CG Air Deluxe (Includes a (1) 3D Mesh Covers) - 22" to 24" Wide	\$910
<input type="checkbox"/>	CGADLX22-24-IC	CG Air Deluxe (Includes an (1) Incontinence Covers) - 22" to 24" Wide	\$910

26" + WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX26+-SC	CG Air Deluxe (Includes a (1) Standard 2-way Stretch Covers) - 26" Wide or more	\$990
<input type="checkbox"/>	CGADLX26+-3D	CG Air Deluxe (Includes a (1) 3D Mesh Covers) - 26" Wide or more	\$1,010
<input type="checkbox"/>	CGADLX26+-IC	CG Air Deluxe (Includes an (1) Incontinence Covers) - 26" Wide or more	\$1,010

REPLACEMENT COVERS

<input type="checkbox"/>	CGASCxxx-x	Standard 2-way Stretch Cover. Specify Dimensions: W _____ " D _____ " H _____ "	\$96
<input type="checkbox"/>	CGA3Dxxx-x	3D Mesh Cover. Specify Dimensions: W _____ " D _____ " H _____ "	\$144
<input type="checkbox"/>	CGAICxxx-x	Incontinence Cover. Specify Dimensions: W _____ " D _____ " H _____ "	\$144

CG Air DELUXE - CUSTOM ORDER MATRIX

WIDTH - FRONT OF CUSHION		
	2 4 6 8 10 12 14 16 18 20 22 24 26 28 30	
DEPTH	2	
	4	
	6	
	8	
	10	
	12	
	14	
	16	
	18	
	20	
	22	
	24	
	26	
	28	
	30	

STEPS TO FOLLOW TO FILL OUT THE MATRIX

1. Specify the requested size of the cushion W _____ " x D _____ "
2. Specify the cell height _____. If different heights are needed, indicate 2, 3 or 4 in each square to show each cell height.
3. Specify the number of valves as well as the valve position by writing a "V". If more than one valve is required, draw a line to show the compartment area.
4. If cells are not required, indicate an "X" in the square.



THE COLORS ARE USED TO SHOW THE DIFFERENT CELL HEIGHT POSSIBILITIES

The CG Air DELUXE can be personalized to the needs of the user as a custom request. The cells are available in heights of 2", 3" and 4". The CG Air DELUXE is a cushion that shapes to the user's body offering optimal pressure distribution and reducing shear forces as well as the channeling of the legs which will increase lateral stability.



CG Air is an Amylor Inc Brand

SIZING CHART

CG Air			OTHER MANUFACTURERS		
WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS	WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS
10"	5	10"	10"	5	8.75"
11"	5	10"	11"	6	10.50"
12"	6	12"	12"	6	10.50"
13"	6	12"	13"	7	12.25"
14"	7	14"	14"	8	14"
15"	7	14"	15"	8	14"
16"	8	16"	16"	9	15.75"
17"	8	16"	17"	9	15.75"
18"	9	18"	18"	10	17.5"
19"	9	18"	19"	10	17.5"
20"	10	20"	20"	11	19.25"
21"	10	20"	21"	12	21"
22"	11	22"	22"	12	21"
23"	11	22"	23"	13	22.75"
24"	12	24"	24"	13	22.75"
25"	12	24"	25"	14	24.5"



The **CG Air** is made of medical grade PVC which offers the following benefits:

- ✓ Reduces pressure on the skin and soft tissues.
- ✓ More stable
- ✓ Better air retention
- ✓ More durable
- ✓ Easier to clean - no retention of smell
- ✓ Meets flammability standards Cal TB 117-2013
- ✓ Best warranty - if defective, cushion replacement for 24 months. See the limited