



▼ ACCOUNT INFORMATION

P.O. # or Quote: _____

Dealer Contact Name: _____

Email Address: _____

Phone Number: _____ Fax Number: _____

Bill to (Name): _____

Bill to (Location): _____

Ship to (Name): _____

Address: _____

City: _____

State: _____ Zip Code: _____

User Name or Number: _____

User Weight: _____

User Height: _____

Seat Cushion Thickness: _____

Knee to Heel Length: Left: _____ Right: _____

The order cannot be processed without the following:

- A user name or number, the client's measurements and weight - Amysystems is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

▼ POWER BASE SELECTION & ELECTONICS

▲ Seat to floor height (STFH) is measured at the front of the seat pan with a 0° Pre-Tilt. Lowest STFH (within ¼") is shown below.

Quickie	▼STFH w/ Tilt	▼STFH w/ Combo
<input type="checkbox"/> Xperience <input type="checkbox"/> 1 <input type="checkbox"/> 2	16 ¼"	18 ¼"
<input type="checkbox"/> Xplore <input type="checkbox"/> 1 <input type="checkbox"/> 2	16 ½"	18 ½"
<input type="checkbox"/> Pulse 6	16 ¾"	18 ¾"
<input type="checkbox"/> QM710 / QM720	TBA	TBA
<input type="checkbox"/> P220 / P222 <input type="checkbox"/> Std. <input type="checkbox"/> Wide	18 ½"	N/A
<input type="checkbox"/> S636 / S646 <input type="checkbox"/> Std <input type="checkbox"/> Wide	17 ¼"	19 ½"

▲ Additional charges may apply if base is equipped with anti-tip cables

▼ Electronics
<input type="checkbox"/> VR2
<input type="checkbox"/> 6 key VR2
<input type="checkbox"/> R-Net / R-Net EL

► Pride Quantum

▲ When a chair is equipped with an 11" seat lift and Q-Logic electronics, reduced speed requires power seating options through drive controls (SWTJSS, SWTJSD or SWTJSM)

▲ Q-Logic NE is not compatible with through joystick activation of the power seating options and does not offer reduce speed with seat elevator.

	▼STFH w/ Tilt	▼STFH w/ Combo
<input type="checkbox"/> Q6 Edge, Edge 2.0	17"	19"
<input type="checkbox"/> J6	17"	19"
<input type="checkbox"/> 4Front	TBA	TBA
<input type="checkbox"/> Rival	TBA	TBA

▼ Electronics	
<input type="checkbox"/> VR2	<input type="checkbox"/> Q-Logic 3
<input type="checkbox"/> 6 key VR2	<input type="checkbox"/> Q-Logic NE
<input type="checkbox"/> Dynamic Shark	<input type="checkbox"/> Q-Logic NE +
<input type="checkbox"/> Q-Logic 1 or 2	

Invacare	▼STFH w/ Tilt	▼STFH w/ Combo
<input type="checkbox"/> FDX	TBA	TBA
<input type="checkbox"/> TDX SP <input type="checkbox"/> GP 22 <input type="checkbox"/> GP 24	16 ½"	19"
<input type="checkbox"/> TDX SR	17"	19 ½"
<input type="checkbox"/> TDX SC (Spree)	15 ½"	17 ½"
<input type="checkbox"/> Storm 3G <input type="checkbox"/> Short <input type="checkbox"/> Long	17"	19 ½"

▼ Electronics	▼ Joysticks
<input type="checkbox"/> Mk690	<input type="checkbox"/> SPJ
<input type="checkbox"/> Mk690ACC	<input type="checkbox"/> MPJ / DISPLAY
<input type="checkbox"/> Mk6TT	
<input type="checkbox"/> Mk6 GTRAC	
<input type="checkbox"/> Linx	

Other Manufacturers	▼STFH w/ Tilt	▼STFH w/ Combo
<input type="checkbox"/> Golden Compass	16 ½"	19"
<input type="checkbox"/> Rovi x3	TBA	TBA
<input type="checkbox"/> IIM - Frontier V6	18 ½"	18 ½"
<input type="checkbox"/> IIM - X8 Extreme	19"	20"
<input type="checkbox"/> Trackmaster MK1	19 ½"	20 ½"
<input type="checkbox"/> Other Power Wheelchair	Specify: _____	

▼ Electronics
<input type="checkbox"/> Golden - Dynamic Shark
<input type="checkbox"/> Rovi - R-Net
<input type="checkbox"/> IIM - PG R-Net
<input type="checkbox"/> Trackmaster - R-net

▼ **INSTALLATION**

<input type="checkbox"/> NOTE-IN-HOUSE	Factory Installation ▲ The Wheelchair will drop-ship from the manufacturer	STD
---	---	-----

▼ **INTERFACE POWER MODULES**

<input type="checkbox"/> PTS50CG	50° CG Power Tilt Interface	\$ 3,495.00
<input type="checkbox"/> PES-C4200	11" ADD-ON Power Elevating Seat Interface Module for transfer assist	\$ 3,495.00

▼ **ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS**

<input type="checkbox"/> SAC	Single (1) Power Seat Function Controller to Interface through R-net Electronics	STD
<input type="checkbox"/> DAC	Dual (2) Power Seat Function Controller to Interface through R-net Electronics	\$ 495.00
With the SAC or DAC, please select one activation and position from the list below		
<input type="checkbox"/> SWTOGGLE	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC) or the Dual (2 Actuators) Power Seat Function Controller (DAC)	<input type="checkbox"/> User <input type="checkbox"/> Attendant STD
<input type="checkbox"/> SWPUSH1	Single Push Button Switch (Single Mode: SAC) ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC)	<input type="checkbox"/> User <input type="checkbox"/> Attendant STD
<input type="checkbox"/> SWPUSH2	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC) or the Dual (2 Actuators) Power Seat Function Controller (DAC)	<input type="checkbox"/> User <input type="checkbox"/> Attendant STD
<input type="checkbox"/> SWTJSS	Single actuator activation through R-net Joystick or Specialty Controls ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC)	\$ 995.00
<input type="checkbox"/> SWTJSD	Dual actuator activation through R-net Joystick or Specialty Controls ▲ Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)	\$ 995.00
<input type="checkbox"/> SWTVR2	Tilt only activation through 6-key PG Drives VR2 Joysticks ▲ Requires the P&G 6-Key VR2 Joystick from Amysystems, Quickie or Pride.	N/C
<input type="checkbox"/> SWTMK6-ACCT	Tilt only activation through Invacare: ▲ Requires the SPJ-ACC, SPJ-PSS or MPJ Joystick with the MK690-ACC controller from Invacare.	<input type="checkbox"/> SPJ-PSS <input type="checkbox"/> SPJ-ACC <input type="checkbox"/> MPJ N/C

▼ **JOYSTICK MOUNTING LOCATION**

<input type="checkbox"/> HML	Install Joystick on Left hand side
<input type="checkbox"/> HMR	Install Joystick on Right hand side

▼ **SEAT SELECTION**

<input type="checkbox"/> Quickie Rehab Seat	<input type="checkbox"/> IIM Rehab Seat	<input type="checkbox"/> Invacare ADJASBA Adult
<input type="checkbox"/> Quickie Sedeo	<input type="checkbox"/> IIM MPS Captain seat	<input type="checkbox"/> Pride Synergy
<input type="checkbox"/> Quickie Eclipse	<input type="checkbox"/> Invacare ADJASBA Pediatric	<input type="checkbox"/> Pride TB2 or TB3
Specify Seat Width: _____	Specify Seat Depth: _____	

▼ **FRONT RIGGINGS**

<input type="checkbox"/> SOMFR	Using the OEM footrests from:	<input type="checkbox"/> Invacare <input type="checkbox"/> Pride	\$ 295.00
<input type="checkbox"/> SOAAF90AO	90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°) ▲ Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB). With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate. Select Foot Platform Range: <input type="checkbox"/> SOFBA7 (7 ½" - 12 ½") <input type="checkbox"/> SOFBA12 (12 ½" - 16 ½") <input type="checkbox"/> SOFBA16 (16 ½" - 19 ½")	<input type="checkbox"/> Innovation in Motion <input type="checkbox"/> Quickie	\$ 225.00
<input type="checkbox"/> SOFUFP	Flip-up one piece Foot Platform		N/C
<input type="checkbox"/> SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each		\$ 63.00
<input type="checkbox"/> SOACPCM	Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each		\$ 63.00

▼ **SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)**

<input type="checkbox"/> SOMHSS	Specific seat to floor height requirement other than standard. Specify at: _____ " (± ¼" Range)	RFQ
--	--	-----

▼ **NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS**

