



**▼ ACCOUNT INFORMATION**

P.O. # or Quote: \_\_\_\_\_

Dealer Contact Name: \_\_\_\_\_

Email Address: \_\_\_\_\_

Phone Number: \_\_\_\_\_ Fax Number: \_\_\_\_\_

Bill to (Name): \_\_\_\_\_

Bill to (Location): \_\_\_\_\_

Ship to (Name): \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_

User Name or Number: \_\_\_\_\_

User Weight: \_\_\_\_\_

User Height: \_\_\_\_\_

Seat Cushion Thickness: \_\_\_\_\_

Knee to Heel Length: Left: \_\_\_\_\_ Right: \_\_\_\_\_

The order cannot be processed without the following:

- A user name or number, the client's measurements and weight - Amysystems is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**▼ POWER BASE SELECTION & ELECTONICS**

▲ Seat to floor height (STFH) is measured at the front of the seat pan with a 0° Pre-Tilt. Lowest STFH (within ¼") is shown below.

Quickie	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics	
<input type="checkbox"/> Xperience <input type="checkbox"/> 1 <input type="checkbox"/> 2	16 ¼"	18 ¼"	<input type="checkbox"/> VR2	
<input type="checkbox"/> Xplore <input type="checkbox"/> 1 <input type="checkbox"/> 2	16 ½"	18 ½"	<input type="checkbox"/> 6 key VR2	
<input type="checkbox"/> Pulse 6	16 ¾"	18 ¾"	<input type="checkbox"/> R-Net / R-Net EL	
<input type="checkbox"/> QM710 / QM720	TBA	TBA		
<input type="checkbox"/> P220 / P222 <input type="checkbox"/> Std. <input type="checkbox"/> Wide	18 ½"	N/A		
<input type="checkbox"/> S636 / S646 <input type="checkbox"/> Std <input type="checkbox"/> Wide	17 ¼"	19 ½"		
▲ Additional charges may apply if base is equipped with anti-tip cables				
Pride Quantum	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics	
▲ When a chair is equipped with an 11" seat lift and Q-Logic electronics, reduced speed requires power seating options through drive controls (SWTJSS, SWTJSD or SWTJSM)				
▲ Q-Logic NE is not compatible with through joystick activation of the power seating options and does not offer reduce speed with seat elevator.				
<input type="checkbox"/> Q6 Edge, Edge 2.0	17"	19"	<input type="checkbox"/> VR2 <input type="checkbox"/> Q-Logic 3	
<input type="checkbox"/> J6	17"	19"	<input type="checkbox"/> 6 key VR2 <input type="checkbox"/> Q-Logic NE	
<input type="checkbox"/> 4Front	TBA	TBA	<input type="checkbox"/> Dynamic Shark <input type="checkbox"/> Q-Logic NE +	
<input type="checkbox"/> Rival	TBA	TBA	<input type="checkbox"/> Q-Logic 1 or 2	
Invacare	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics	▼ Joysticks
<input type="checkbox"/> FDX	TBA	TBA	<input type="checkbox"/> Mk690	<input type="checkbox"/> SPJ
<input type="checkbox"/> TDX SP <input type="checkbox"/> GP 22 <input type="checkbox"/> GP 24	16 ½"	19"	<input type="checkbox"/> Mk690ACC	<input type="checkbox"/> MPJ / DISPLAY
<input type="checkbox"/> TDX SR	17"	19 ½"	<input type="checkbox"/> Mk6TT	
<input type="checkbox"/> TDX SC (Spree)	15 ½"	17 ½"	<input type="checkbox"/> Mk6 GTRAC	
<input type="checkbox"/> Storm 3G <input type="checkbox"/> Short <input type="checkbox"/> Long	17"	19 ½"	<input type="checkbox"/> Linx	
Other Manufacturers	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics	
<input type="checkbox"/> Golden Compass	16 ½"	19"	<input type="checkbox"/> Golden - Dynamic Shark	
<input type="checkbox"/> Rovi x3	TBA	TBA	<input type="checkbox"/> Rovi - R-Net	
<input type="checkbox"/> IIM - Frontier V6	18 ½"	18 ½"	<input type="checkbox"/> IIM - PG R-Net	
<input type="checkbox"/> IIM - X8 Extreme	19"	20"		
<input type="checkbox"/> Trackmaster MK1	19 ½"	20 ½"	<input type="checkbox"/> Trackmaster - R-net	
<input type="checkbox"/> Other Power Wheelchair	Specify: _____			

## ▼ INSTALLATION

<input type="checkbox"/> NOTE-IN-HOUSE	Factory Installation ▲ The Wheelchair will drop-ship from the manufacturer	STD
--	---	-----

## ▼ POWER OPTIONS & SEATING

<input type="checkbox"/> SAABA	Amyseat Rehab Seat - Seat Frame only	\$ 1,495.00
<input type="checkbox"/> SAPTN	CG Power Tilt ▲ This option includes the Amyseat Rehab Seat (SAABA)	\$ 6,495.00
<input type="checkbox"/> SAPRN	Power Recline with Integrated Shear Reduction (ISR) ▲ This option includes the Amyseat Rehab Seat (SAABA). It is recommended to use the NXT Back Support option with the power or manual reclines (NXTFHTAC)	\$ 7,995.00
<input type="checkbox"/> SAPTRN	CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR) ▲ This option includes the Amyseat Rehab Seat (SAABA). It is recommended to use the NXT Back Support option with the power or manual reclines (NXTFHTAC)	\$ 11,495.00
<input type="checkbox"/> SAPES	11" Power Elevating Seat Module for transfer assist ▲ This option includes the Amyseat Rehab Seat (SAABA)	\$ 4,495.00
<input type="checkbox"/> ADDPES	<b>Module - Add an 11" Power Elevating Seat to either SAPTX, SAPRX or SAPTRX</b>	\$ 2,995.00
<input type="checkbox"/> ADDMR	<b>Module - Manual Recline with Integrated Shear Reduction (ISR)</b> ▲ It is recommended to use the NXT Back Support option with the power or manual reclines (NXTFHTAC). The maximum user weight capacity with this option is 250 lbs	\$ 1,495.00

## ▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/> SAC	Single (1) Power Seat Function Controller to Interface through R-net Electronics	STD
<input type="checkbox"/> DAC	Dual (2) Power Seat Function Controller to Interface through R-net Electronics	\$ 495.00
<b>With the SAC or DAC, please select one activation and position from the list below</b>		
<input type="checkbox"/> SWTOGGLE	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User <input type="checkbox"/> Attendant	STD
▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC) or the Dual (2 Actuators) Power Seat Function Controller (DAC)		
<input type="checkbox"/> SWPUSH1	Single Push Button Switch (Single Mode: SAC) <input type="checkbox"/> User <input type="checkbox"/> Attendant	STD
▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC)		
<input type="checkbox"/> SWPUSH2	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) <input type="checkbox"/> User <input type="checkbox"/> Attendant	STD
▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC) or the Dual (2 Actuators) Power Seat Function Controller (DAC)		
<input type="checkbox"/> SWTJSS	Single actuator activation through R-net Joystick or Specialty Controls	\$ 995.00
▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC)		
<input type="checkbox"/> SWTJSD	Dual actuator activation through R-net Joystick or Specialty Controls	\$ 995.00
▲ Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)		
<input type="checkbox"/> SWTVR2	Tilt only activation through 6-key PG Drives VR2 Joysticks	N/C
▲ Requires the P&G 6-Key VR2 Joystick from Amysystems, Quickie or Pride.		
<input type="checkbox"/> SWTMK6-ACCT	Tilt only activation through Invacare: <input type="checkbox"/> SPJ-PSS <input type="checkbox"/> SPJ-ACC <input type="checkbox"/> MPJ	N/C
▲ Requires the SPJ-ACC, SPJ-PSS or MPJ Joystick with the MK690-ACC controller from Invacare.		

## ▼ ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/> MAC	Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	\$ 2,195.00
<b>With the MAC, please select one activation and position from the list below</b>		
<input type="checkbox"/> SWMULTI	4-Way Toggle Switch (Up to 4 Activations) <input type="checkbox"/> User <input type="checkbox"/> Attendant	N/C
▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)		
<input type="checkbox"/> SWMULTIT	Triple Toggle Switch (Up to 5 Activations) <input type="checkbox"/> User <input type="checkbox"/> Attendant	N/C
▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)		
<input type="checkbox"/> SWPUSH4	4-Way Push Button Switch <input type="checkbox"/> User <input type="checkbox"/> Attendant	N/C
▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)		
<input type="checkbox"/> SWPUSH5	5-Way Push Button Switch <input type="checkbox"/> User <input type="checkbox"/> Attendant	N/C
▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)		
<input type="checkbox"/> SWTJSM	Multiple Actuator activation through R-net Joystick or Specialty Controls	STD
▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)		

## ▼ SWITCH ACTIVATION HAND MOUNT

<input type="checkbox"/> TSML	Install Switch on LH side	N/C
<input type="checkbox"/> TSMR	Install Switch on RH side	N/C

## ▼ JOYSTICK MOUNT

<input type="checkbox"/> SOSWJSM	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis <input type="checkbox"/> Set-up with standard installation to swing-away sideways <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards	\$ 364.00
----------------------------------	--	-----------

## ▼ JOYSTICK MOUNT SUPPORT

<input type="checkbox"/> JSMFH	Joystick mount support Bracket	STD
--------------------------------	--------------------------------	-----

## ▼ JOYSTICK HAND MOUNT

<input type="checkbox"/> HML	Install Joystick on Left hand side
<input type="checkbox"/> HMR	Install Joystick on Right hand side

▼ **PRE-SET SEAT ANGLE (For Non-Reclining Backs)**

**SOBAA** Please select one:  86°  94°  102°  110°  118°  90°  98°  106°  114°  122° N/C

▼ **PEDIATRIC ADJUSTABLE SEAT WIDTH RANGE BETWEEN 11" AND 16"**

**SOWN11-16** Width adjustment from 11" to 16" and Depth range of 11" to 16" N/C  
 ▲ Seat to Floor Height increases by 1/4" N/C

11" (SOWN11)  13" (SOWN13)  15" (SOWN15)  
 12" (SOWN12)  14" (SOWN14)  16" (SOWN16)

▼ **PEDIATRIC ADJUSTABLE SEAT DEPTH RANGE BETWEEN 11" AND 16"**

11" (SODN11)  13" (SODN13)  15" (SODN15)  
 12" (SODN12)  14" (SODN14)  16" (SODN16)

▼ **STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"**

**SOWN15-20-19** Width adjustment from 15" to 20" with a maximum Depth of 19" N/C  
 **SOWN15-20-22** Width adjustment from 15" to 20" with a maximum Depth of 22" N/C

15" (SOWN15)  17" (SOWN17)  19" (SOWN19S)  
 16" (SOWN16)  18" (SOWN18)  20" (SOWN20S)

▼ **LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24"**

**SOWN19-24-22** Width adjustment from 19" to 24" and Depth range of 18" to 22" \$ 430.00

19" (SOWN19L)  21" (SOWN21)  23" (SOWN23)  
 20" (SOWN20L)  22" (SOWN22)  24" (SOWN24)

▼ **CUSTOM WIDTH ABOVE 25"**

**SOWCUS** Custom Seat Width Specify: \_\_\_\_\_ \$ 465.00

▼ **ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22"**

15" (SODN15)  18" (SODN18)  21" (SODN21)  
 16" (SODN16)  19" (SODN19)  22" (SODN22)  
 17" (SODN17)  20" (SODN20)

▼ **CUSTOM DEPTH ABOVE 22"**

**SODCUS** Custom Seat Depth Specify: \_\_\_\_\_ \$ 250.00

▼ **BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)**

16" (SOBCH16)  22" (SOBCH22) 22" is the minimum back cane height required with a power or manual recline.  
 18" (SOBCH18)  24" (SOBCH24)  
 20" (SOBCH20)  26" (SOBCH26)

▼ **CUSTOM BACK CANE HEIGHT ABOVE 27"**

**SOBHCUS** Custom Back Cane Height Specify: \_\_\_\_\_ \$ 250.00

▼ **BACK SUPPORT**

**STDSB** Standard contour back support 15" to 24" Wide with a maximum shell height of 22" \$ 550.00  
 ▲ The minimum space between the back and seat pan is 3" with a fixed back and 6" with a reclining back

**NXTFHTAC** NXT Optima FHT AC Fixed Thoracic Back Support (15" to 21" Wide) \$ 715.00

**NXTFHTAC** NXT Optima FHT AC Fixed Thoracic Back Support (22" to 24" Wide) \$ 762.00

▲ The actual back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 3" with a fixed back and 6" with a reclining back

**STDSBO** Omission of the Contour Back Support Specify Model of Back Support: \_\_\_\_\_ N/C

▲ It is recommended to use the Back Support option with the power or manual recline.

▼ **BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

16" (SOBH16)  24" (SOBH24)  
 18" (SOBH18)  26" (SOBH26)  
 20" (SOBH20)  28" (SOBH28) This height is not available with the NXT Optima option (NXTFHTAC)  
 22" (SOBH22) 22" is the minimum back height required with a power or manual recline.

▼ SEAT CUSHIONS

► CG Air - THERAPEUTIC AIR CUSHION

<input type="checkbox"/> Single Valve 15" to 21" Wide		CELL HEIGHT <input type="checkbox"/> 2"	
<input type="checkbox"/> Single Valve 22" to 24" Wide		<input type="checkbox"/> 3"	
<input type="checkbox"/> Custom Cushion (Please use the CG Air Order Form)		<input type="checkbox"/> 4"	
<input type="checkbox"/> CGA1VSCxxxx-x	CG Air Single Valve (with Standard 2-way Stretch Cover) - NEW		\$ 630.00
<input type="checkbox"/> CGA1V3Dxxxx-x	CG Air Single Valve (with 3D Mesh Cover) - NEW		\$ 740.00
<input type="checkbox"/> CGA1VICxxxx-x	CG Air Single Valve (with Incontinence Cover) - NEW		\$ 740.00
▲ The actual cushion sizes only come in even numbers for width and depth measurements.			

► NXT BIO-FIT CUSHION

<input type="checkbox"/> NXTBIOFIT	NXT BIO FIT Cushion (15" to 21" Wide)		\$ 455.00
<input type="checkbox"/> NXTBIOFIT	NXT BIO FIT Cushion (22" to 24" Wide)		\$ 482.00
▲ The actual cushion sizes only come in even numbers for width and depth measurements.			

▼ CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)

<input type="checkbox"/> SOAFBFLH	Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments)		STD
<input type="checkbox"/> SOAFBFRH	Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments)		STD
<input type="checkbox"/> SOAFBDLH	Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments)		N/C
<input type="checkbox"/> SOAFBDRH	Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments)		N/C

▼ DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)

<input type="checkbox"/> SOADPPLH	Pediatric Dual Post Height Adjustable - Left Side (Height Adjustable from 7 ¼" - 10 ¼")		\$ 175.00
<input type="checkbox"/> SOADPPRH	Pediatric Dual Post Height Adjustable - Right Side (Height Adjustable from 7 ¼" - 10 ¼")		\$ 175.00
<input type="checkbox"/> SOADPFLH	Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")		\$ 175.00
<input type="checkbox"/> SOADPFRH	Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")		\$ 175.00
<input type="checkbox"/> SOADPDLH	Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")		\$ 175.00
<input type="checkbox"/> SOADPDRH	Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 10 ¼" - 15 ¼")		\$ 175.00

▼ RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)

<input type="checkbox"/> SOANIFLH	Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)		N/C
<input type="checkbox"/> SOANIFRH	Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments)		N/C
<input type="checkbox"/> SOANIDLH	Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)		N/C
<input type="checkbox"/> SOANIDRH	Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments)		N/C

▼ ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: \_\_\_\_\_ "

▲ The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post

▼ ARMPADS

► Foam Type Selection

<input type="checkbox"/> SOFAPD	Desk Length Standard (10")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	STD
<input type="checkbox"/> SOFAPF	Full Length Standard (14")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	STD
<input type="checkbox"/> SOWAPD	Desk Length Waterfall (10")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$ 30.00
▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports			
<input type="checkbox"/> SOWAPF	Full Length Waterfall (14")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$ 30.00
▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports			

► Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/> SOGELD	Desk Length Hybrid Gel and Foam Flat armpad (10")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$ 60.00
<input type="checkbox"/> SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$ 60.00
<input type="checkbox"/> SOWGELD	Desk Length Hybrid Gel and Foam Waterfall armpad (10")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$ 95.00
▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports			
<input type="checkbox"/> SOWGELF	Full Length Hybrid Gel and Foam Waterfall armpad (14")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$ 95.00
▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports			

► Specialty Type Selection

<input type="checkbox"/> SOCAR	Arm Trough with Hand Support (18" x 4 ½")	Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$ 195.00
▲ Arm Trough Armpads may interfere with Armrests Flipping back.			

## ▼ FOOTRESTS

<input type="checkbox"/>	<b>SOPAAF</b>	Pediatric Swing-Away Footrests ▲ Depending on seat configuration, 70°, 80° and 90° might not be available. Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70° <input type="checkbox"/> 80° <input type="checkbox"/> 90° Select Footrest Range: <input type="checkbox"/> SOFLP5 (5 ½" - 7 ½") <input type="checkbox"/> SOFLP8 (8 ½" - 10 ½") <input type="checkbox"/> SOFLP10 (10 ½" - 15 ½")	\$	295.00
<input type="checkbox"/>	<b>SOSA AF</b>	Adult Swing-Away Footrests ▲ This option is not available with the pediatric seat frame. Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70° Select Footrest Range: <input type="checkbox"/> SOFLA10 (10 ½" - 15 ½") <input type="checkbox"/> SOFLA14 (14 ½" - 19 ½")	\$	245.00
<input type="checkbox"/>	<b>SOAAF</b>	Adult <b>Heavy Duty</b> Swing-Away footrests ▲ This option is not available with the pediatric seat frame. Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70° Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 ½" - 15 ½") <input type="checkbox"/> SOFLHD14 (14 ½" - 19 ½")	\$	345.00
<input type="checkbox"/>	<b>ADDPL</b>	Adult Power Elevating & Articulating Legrests ▲ This option is not available with the pediatric seat frame. This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function. Please Select Function: <input type="checkbox"/> Individual <input type="checkbox"/> Combined Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 ½" - 15 ½") <input type="checkbox"/> SOFLPW14 (14 ½" - 19 ½")	\$	2,795.00

## ▼ FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	<b>SOSFAF</b>	Flip-up Footplates		STD
<input type="checkbox"/>	<b>SOAFP</b>	Angle Adjustable & Flip-up Footplates - Each	\$	62.00
<input type="checkbox"/>	<b>SOAHL</b>	Heel Loops - Each	\$	20.00
<input type="checkbox"/>	<b>SOACPFR</b>	Calf Pads for Adult Heavy-Duty Swing-Away footrests ( <b>SOAAF</b> ) with adjustable hardware - Each	\$	63.00

## ▼ FOOTREST LENGTH ADJUSTMENT AT FACTORY

Please set the Footrest Length at: \_\_\_\_\_ "

The Footrest Length is taken from the top of the seat pan to the top of the footplate.

## ▼ CENTER MOUNTS

<input type="checkbox"/>	<b>SOAAF90</b>	90° Fixed Center Mount Interface Bracket ( <b>Offers Manual adjustability from 95° to 70°</b> ) ▲ Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB). With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate. Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 ½" - 16 ½") <input type="checkbox"/> SOFBA16 (16 ½" - 19 ½")	\$	225.00
<input type="checkbox"/>	<b>ADDPF</b>	90° Power Center Mount Elevating & Articulating Foot Platform ▲ This option includes calf pads. The power center mount counts as one power function. Please Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB). This option is not available with the 12 ½" Select Foot Platform Range: <input type="checkbox"/> SOFBPW12 (12 ½" - 16 ½") <input type="checkbox"/> SOFBPW16 (16 ½" - 20 ½")	\$	2,795.00

## ▼ FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	<b>SOFUFP</b>	Flip-up one piece Foot Platform		N/C
<input type="checkbox"/>	<b>SOAFPFB</b>	Angle Adjustable & Flip-up Split Footplates for center mount - Each	\$	63.00
<input type="checkbox"/>	<b>SOACPCM</b>	Calf Pads for 90° Fixed Center Mount ( <b>SOAAF90</b> ) with adjustable hardware - Each	\$	63.00

## ▼ FOOT PLATFORM LENGTH ADJUSTMENT AT FACTORY

Please set the Foot Platform Length at: \_\_\_\_\_ "

The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

## ▼ SEAT POSITIONING

<input type="checkbox"/>	<b>SOAUTO</b>	55" Seat Positioning Strap		STD
<input type="checkbox"/>	<b>SOAUTO80</b>	80" Seat Positioning Strap	\$	75.00
<input type="checkbox"/>	<b>SOBPHP2C</b>	Bodypoint padded Hip Belt (Two Point - Center Pull)	\$	164.00
<input type="checkbox"/>	<b>SOBPCS</b>	Bodypoint padded Chest Strap	\$	235.00

## ▼ HEADREST

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

<input type="checkbox"/>	<b>SONUAHR</b>	Using headrest mounting hardware other than SOMAH located below. Specify: _____		N/C
<input type="checkbox"/>	<b>SOFHRS</b>	Small Curved Headrest Pad 4" x 9"	\$	200.00
<input type="checkbox"/>	<b>SOFHRM</b>	Medium Curved Headrest Pad 5" x 10"	\$	200.00
<input type="checkbox"/>	<b>SOFHRL</b>	Large Curved Headrest Pad 6" x 11"	\$	200.00
<input type="checkbox"/>	<b>SOMAH</b>	Multi-axis Removable Mounting Hardware	\$	160.00
<input type="checkbox"/>	<b>SOHRBC</b>	Bracket to Mount on the Contour Back Support ( <b>Without Recline</b> )	\$	100.00
<input type="checkbox"/>	<b>SOHRPR</b>	Bracket to Mount on the Seat Frame ( <b>For Reclining Backs</b> )	\$	100.00

### ▼ LATERAL TRUNK SUPPORTS

Please select one of the lateral pads and the adjustable & removable hardware.

<input type="checkbox"/>	<b>SOSLS3</b>	Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"	\$	100.00
<input type="checkbox"/>	<b>SOSLM4</b>	Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6"	\$	100.00
<input type="checkbox"/>	<b>SOSLL5</b>	Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"	\$	100.00
<input type="checkbox"/>	<b>SOSLHM</b>	Adjustable Swing-Away Hardware Mounts - Each	\$	200.00

### ▼ HIP & THIGH SUPPORTS

Please select one of the hip & thigh support pads and the adjustable & removable hardware.

<input type="checkbox"/>	<b>SOHTSS10</b>	Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"	\$	115.00
<input type="checkbox"/>	<b>SOHTSM12</b>	Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5"	\$	115.00
<input type="checkbox"/>	<b>SOHTSL14</b>	Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"	\$	115.00
<input type="checkbox"/>	<b>SOHTSARH</b>	Adjustable and Removable Hip & Thigh Support Hardware - Each	\$	215.00

### ▼ ACCESSORIES

<input type="checkbox"/>	<b>SOCH</b>	Cup Holder (Self-Leveling)	\$	85.00
<input type="checkbox"/>	<b>SOCPH</b>	Ram Mounts X-Grip Cell Phone Holder	\$	110.00
<input type="checkbox"/>	<b>SOTHS</b>	Ram Mounts X-Grip Small Tablet Holder	\$	160.00
<input type="checkbox"/>	<b>SOTHL</b>	Ram Mounts X-Grip Large Tablet Holder with Swivel	\$	250.00

### ▼ BACK PACK AND ARMREST SIDE POUCHES

<input type="checkbox"/>	<b>SOBP</b>	Back Pack (Comes with durable straps that can be hooked around the back canes)	\$	88.00
<input type="checkbox"/>	<b>SOAP</b>	Armrest Side Pouch (Comes with durable straps that can be hooked around the armrest tubes) - Order <input type="checkbox"/> 1 <input type="checkbox"/> 2	\$	62.00

### ▼ SWITCHES AND OPTIONS

<input type="checkbox"/>	<b>SSBS</b>	Buddy Button Switch	\$	90.00
<input type="checkbox"/>	<b>SSES</b>	Egg Switch-Black	\$	90.00
<input type="checkbox"/>	<b>SSML</b>	Ultra Light Switch	\$	145.00

### ▼ CUSTOM SECTION

Lead times will increase depending on the custom order pipeline at the time the purchase order is received. Please call for Quote and specifications for these options.

#### ▼ O<sup>2</sup> Holder

<input type="checkbox"/>	<b>SUPOXY</b>	Custom Non-Articulating Oxygen Tank Holder ▲ Not available with a power tilt and/or power recline	\$	295.00
<input type="checkbox"/>	<b>SUPOXYA</b>	Custom Articulating Oxygen Tank Holder ▲ Required with a power tilt and/or power recline	\$	1,195.00

- D-Size Oxygen Tank (\* 22" Tall x 4.3" Diameter) \* The height includes the valve
- E-Size Oxygen Tank (\* 31" Tall x 4.3" Diameter) \* The height includes the valve
- MLH (Mount on Left)
- MRH (Mount on Right)

#### ▼ VENT TRAY (Please complete in Full)

Amysystems does not supply the battery, charger or the wiring for the Vent. These items are normally provided by the Vent Manufacturer.

<input type="checkbox"/>	<b>SOVENT</b>	Custom Vent Tray		
		<input type="checkbox"/> <b>Articulating Vent Box</b> ▲ Required with Power Tilt and/or Recline	\$	1,360.00
		<input type="checkbox"/> <b>Fixed Vent Box</b>	\$	400.00
		Specify Vent Model: _____		
<input type="checkbox"/>	<b>SOVTBE</b>	Custom Tray for Ventilator Battery	\$	495.00
		<input type="checkbox"/> The Vent is being powered by 1 U1 Battery		
		<input type="checkbox"/> The Vent is being powered by 1 22NF Battery		
		<input type="checkbox"/> The Vent is being powered by 1 Group 24 Battery		
		<input type="checkbox"/> The Vent is being powered by 1 Group 34 Battery		

▼ **CUSTOM SECTION - Continued**

Lead times will increase depending on the custom order pipeline at the time the purchase order is received. Please call for Quote and specifications for these options.

▼ **OTHER CUSTOM OPTIONS (Request For Quote)**

Please specify details and specifications for these options.

<input type="checkbox"/> <b>PAT</b>	Power Anterior Tilt	RFQ
<input type="checkbox"/> <b>CUSTOM</b>	4-Point wheelchair tie downs ▲ 4-Point wheelchair transport brackets are not crash-tested in accordance with WC-19, Use these transport brackets only to secure an unoccupied wheelchair during transport.	RFQ
<input type="checkbox"/> <b>CUSTOM</b>	Power Flip-up Foot Platform (Please call CSR for availability with configuration and user weight capacity)	RFQ
<input type="checkbox"/> <b>CUSTOM</b>	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan	RFQ
<input type="checkbox"/> <b>CUSTOM</b>	Endomorphic Back Canes (Pair) <input type="checkbox"/> 1" per side towards the inside = 2" in <input type="checkbox"/> 2" per side towards the inside = 4" in	RFQ
<input type="checkbox"/> <b>CUSTOM</b>	Transfer Handles (Pair)	RFQ
<input type="checkbox"/> <b>CUSTOM # 1</b>	Other custom options. Specify in notes section bellow.	RFQ
<input type="checkbox"/> <b>CUSTOM # 2</b>	Other custom options. Specify in notes section bellow.	RFQ

▼ **SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)**

<input type="checkbox"/> <b>SOMHSS</b>	Specific seat to floor height requirement other than standard. Specify at: _____" ( $\pm$ 1/4" Range)	RFQ
--	--	-----

▼ **NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS**

