

▼ ACCOUNT INFORMATION

P.O. # or Quote: _____
 Dealer Contact Name: _____
 Email Address: _____
 Phone Number: _____ Fax Number: _____
 Bill to (Name): _____
 Bill to (Location): _____
 Ship to (Name): _____
 Address: _____
 City: _____
 State/Province: _____ Zip Code: _____
 User Name or Number: _____
 User Weight: _____
 User Height: _____
 Seat Cushion Thickness: _____
 Knee to Heel Length: Left: _____ Right: _____

The order cannot be processed without the following:

- A user name or number, the client's measurements and weight - Amysystems is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES

▼ POWER BASE SELECTION & ELECTONICS

▲ Seat to floor height (STFH) is measured at the front of the seat pan with a 0° Pre-Tilt. Lowest STFH (within ¼") is shown below.

Quickie	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics
<input type="checkbox"/> Xperience <input type="checkbox"/> 1 <input type="checkbox"/> 2	16 ¼"	18 ¼"	<input type="checkbox"/> VR2
<input type="checkbox"/> Xplore <input type="checkbox"/> 1 <input type="checkbox"/> 2	16 ½"	18 ½"	<input type="checkbox"/> 6 key VR2
<input type="checkbox"/> Pulse 6	16 ¾"	18 ¾"	<input type="checkbox"/> R-Net 90 <input type="checkbox"/> R-Net 120
<input type="checkbox"/> QM Series -Specify: _____	TBA	TBA	
<input type="checkbox"/> P220 / P222 <input type="checkbox"/> Std. <input type="checkbox"/> Wide	18 ½"	N/A	
<input type="checkbox"/> S636 / S646 <input type="checkbox"/> Std <input type="checkbox"/> Wide	17 ¼"	19 ½"	

▲ Additional charges may apply if base is equipped with anti-tip cables

► Pride Quantum

▲ When a chair is equipped with a seat elevate with any Q-Logic electronics, reduced speed requires power seating options through drive controls (SWTJSS, SWTJSD or SWTJSM)

▲ Q-Logic NE is not compatible with through joystick activation of the power seating options and does not offer reduced speed with seat elevator. Only Drive-Lock-out is available.

▲ Q-Logic NE+ is only available for a maximum of one actuator for through the Joystick activation.

▲ Dynamic Shark is only available for a maximum of two actuators with external switches. For through the joystick activation, only one actuator can be used and it requires one of the following joysticks be used: Shark PSS or Shark ACC joysticks

	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics
<input type="checkbox"/> Q6 Edge, Edge 2.0	17"	19"	<input type="checkbox"/> VR2 <input type="checkbox"/> Q-Logic 3
<input type="checkbox"/> J6	17"	19"	<input type="checkbox"/> 6 key VR2 <input type="checkbox"/> Q-Logic NE
<input type="checkbox"/> 4Front	TBA	TBA	<input type="checkbox"/> Dynamic Shark <input type="checkbox"/> Q-Logic NE +
<input type="checkbox"/> Rival	TBA	TBA	<input type="checkbox"/> Q-Logic 1 or 2

Invacare	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics	▼ Joysticks
<input type="checkbox"/> FDX	TBA	TBA	<input type="checkbox"/> Mk690	<input type="checkbox"/> SPJ
<input type="checkbox"/> TDX SP <input type="checkbox"/> GP 22 <input type="checkbox"/> GP 24	16 ½"	19"	<input type="checkbox"/> Mk690ACC	<input type="checkbox"/> MPJ / DISPLAY
<input type="checkbox"/> TDX SR	17"	19 ½"	<input type="checkbox"/> Mk6TT	
<input type="checkbox"/> TDX SC (Spree)	15 ½"	17 ½"	<input type="checkbox"/> Mk6 GTRAC	
<input type="checkbox"/> Storm 3G <input type="checkbox"/> Short <input type="checkbox"/> Long	17"	19 ½"		

Other Manufacturers	▼STFH w/ Tilt	▼STFH w/ Combo	▼ Electronics
<input type="checkbox"/> IIM - Frontier V6	18 ½"	18 ½"	<input type="checkbox"/> IIM - PG R-Net
<input type="checkbox"/> IIM - X8 Extreme	19"	20"	
<input type="checkbox"/> Other Power Wheelchair Specify: _____			

▼ INSTALLATION

- NOTE-IN-HOUSE** Factory Installation STD
 ▲ The Wheelchair must be drop-shipped from the manufacturer or Delivered to Amysystems from the dealer

▼ INTERFACE POWER MODULES

- PTS50CG** 50° CG Power Tilt Interface \$3,495
 PES-300 11" ADD-ON Power Elevating Seat Interface Module for transfer assist \$3,495

▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

- SOLO** Single (1) Power Seat Function Controller - Includes an integrated toggle switch (Dual mode) STD
 DUO Dual (2) Power Seat Function Controller - Includes an integrated toggle switch (Single Mode) \$495
With the SOLO or DUO, it is possible to add an additional separate single or dual push button switch.
- SWPUSH1** Single Push Button Switch (Single Mode: SOLO only) **User** **Attendant** STD
 ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SOLO)
- SWPUSH2** Dual Push Button Switch (Single Mode: DUO - Dual Mode: SOLO) **User** **Attendant** STD
 ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SOLO) or the Dual (2 Actuators) Power Seat Function Controller (DUO)
- SWTJSS** Single actuator activation through the Joystick or Specialty Controls \$995
 ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SOLO)
- SWTJSD** Dual actuator activation through the Joystick or Specialty Controls \$995
 ▲ Requires the Dual (2 Actuators) Power Seat Function Controller (DUO)
- SWTVR2** Tilt only activation through 6-key PG Drives VR2 Joysticks N/C
 ▲ Requires the P&G 6-Key VR2 Joystick from Amysystems, Quickie or Pride.
- EPSEL** Tilt only activation through PG Drives R-Net EL Joysticks N/C
 ▲ Requires the P&G R-Net EL Joystick from Amysystems or Quickie
- SWTMK6-ACCT** Tilt only activation through Invacare: **SPJ-PSS** **SPJ-ACC** **MPJ** N/C
 ▲ Requires the SPJ-ACC, SPJ-PSS or MPJ Joystick with the MK690-ACC controller from Invacare.

▼ JOYSTICK MOUNTING LOCATION

- HML** Install Joystick on Left hand side
 HMR Install Joystick on Right hand side

▼ SEAT SELECTION

- | | | |
|---|---|---|
| <input type="checkbox"/> Quickie Rehab Seat | <input type="checkbox"/> IIM Rehab Seat | <input type="checkbox"/> Invacare ADJASBA Adult |
| <input type="checkbox"/> Quickie Sedeo | <input type="checkbox"/> IIM MPS Captain seat | <input type="checkbox"/> Pride Synergy |
| <input type="checkbox"/> Quickie Eclipse | <input type="checkbox"/> Invacare ADJASBA Pediatric | <input type="checkbox"/> Pride TB2 or TB3 |
- Specify Seat Width: _____ Specify Seat Depth: _____

▼ SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)

- SOMHSS** Specific seat to floor height requirement other than standard. RFQ
 Specify at: _____ " (± ¼" Range)

▼ NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS



www.amysystems.com / info@amysystems.com / Phone: +1 450 424-0288 / Fax: +1 450 424-7211