



DECEMBER 2018

**▼ ACCOUNT INFORMATION**

P.O. # or Quote: \_\_\_\_\_  
 Dealer Contact Name: \_\_\_\_\_  
 Email Address: \_\_\_\_\_  
 Phone Number: \_\_\_\_\_ Fax Number: \_\_\_\_\_  
 Bill to (Name): \_\_\_\_\_  
 Bill to (Location): \_\_\_\_\_  
 Ship to (Name): \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_  
 State/Province: \_\_\_\_\_ Zip Code: \_\_\_\_\_  
 User Name or Number: \_\_\_\_\_  
 User Weight: \_\_\_\_\_  
 User Height: \_\_\_\_\_  
 Seat Cushion Thickness: \_\_\_\_\_  
 Knee to Heel Length: Left: \_\_\_\_\_ Right: \_\_\_\_\_

The order cannot be processed without the following:

- A user name or number, the client's measurements and weight
- Amysystems is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amysystems summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES**

**▼ POWER BASE SELECTION & ELECTONICS**

▲ Seat to floor height (STFH) is measured at the front of the seat pan with a 0° Pre-Tilt. Lowest STFH (within ¼") is shown below.

| Quickie  | ▼STFH w/ Tilt | ▼STFH w/ Combo | ▼ Electronics  |
|--|---------------|----------------|--|
| <input type="checkbox"/> Xperience <input type="checkbox"/> 1 <input type="checkbox"/> 2         | 16 ¼"         | 18 ¼"          | <input type="checkbox"/> VR2   |
| <input type="checkbox"/> Xplore <input type="checkbox"/> 1 <input type="checkbox"/> 2            | 16 ½"         | 18 ½"          | <input type="checkbox"/> 6 key VR2                                   |
| <input type="checkbox"/> Pulse 6   | 16 ¾"         | 18 ¾"          | <input type="checkbox"/> R-Net 90 <input type="checkbox"/> R-Net 120 |
| <input type="checkbox"/> QM Series -Specify: _____   | TBA           | TBA            |  |
| <input type="checkbox"/> P220 / P222 <input type="checkbox"/> Std. <input type="checkbox"/> Wide | 18 ½"         | N/A            |  |
| <input type="checkbox"/> S636 / S646 <input type="checkbox"/> Std <input type="checkbox"/> Wide  | 17 ¼"         | 19 ½"          |  |

▲ Additional charges may apply if base is equipped with anti-tip cables

**► Pride Quantum**

▲ When a chair is equipped with a seat elevate with any Q-Logic electronics, reduced speed requires power seating options through drive controls (SWTJSS, SWTJSD or SWTJSM)

▲ Q-Logic NE is not compatible with through joystick activation of the power seating options and does not offer reduced speed with seat elevator. Only Drive-Lock-out is available.

Q-Logic NE+ is only available for a maximum of two actuators for through the Joystick activation. It is only compatible with Tilt / Power Elevating Center Mount or Elevate / Power Elevating Center Mount. Multiple power use is possible with the MAC box but does not offer reduced speed with the seat elevate. Only Drive-Lock-out is available with Seat Elevate, and only external switch activations can be used (Not available with activation through the Joystick).

▲ Dynamic Shark is only available for a maximum of two actuators with external switches. For through the joystick activation, only one actuator can be used and it requires one of the following joysticks be used: Shark PSS or Shark ACC joysticks

|  | ▼STFH w/ Tilt | ▼STFH w/ Combo | ▼ Electronics                           |                                       |
|--|---------------|----------------|---|---------------------------------------|
| <input type="checkbox"/> Q6 Edge, Edge 2.0 | 17"           | 19"            | <input type="checkbox"/> VR2            | <input type="checkbox"/> Q-Logic 3    |
| <input type="checkbox"/> J6                | 17"           | 19"            | <input type="checkbox"/> 6 key VR2      | <input type="checkbox"/> Q-Logic NE   |
| <input type="checkbox"/> 4Front            | TBA           | TBA            | <input type="checkbox"/> Dynamic Shark  | <input type="checkbox"/> Q-Logic NE + |
| <input type="checkbox"/> Rival             | TBA           | TBA            | <input type="checkbox"/> Q-Logic 1 or 2 |                                       |

| Invacare   | ▼STFH w/ Tilt | ▼STFH w/ Combo | ▼ Electronics                      | ▼ Joysticks                            |
|--|---------------|----------------|------------------------------------|--|
| <input type="checkbox"/> FDX   | TBA           | TBA            | <input type="checkbox"/> Mk690     | <input type="checkbox"/> SPJ           |
| <input type="checkbox"/> TDX SP <input type="checkbox"/> GP 22 <input type="checkbox"/> GP 24  | 16 ½"         | 19"            | <input type="checkbox"/> Mk690ACC  | <input type="checkbox"/> MPJ / DISPLAY |
| <input type="checkbox"/> TDX SR  | 17"           | 19 ½"          | <input type="checkbox"/> Mk6TT     |  |
| <input type="checkbox"/> TDX SC (Spree)  | 15 ½"         | 17 ½"          | <input type="checkbox"/> Mk6 GTRAC |  |
| <input type="checkbox"/> Storm 3G <input type="checkbox"/> Short <input type="checkbox"/> Long | 17"           | 19 ½"          |                                    |  |

| Other Manufacturers                             | ▼STFH w/ Tilt  | ▼STFH w/ Combo | ▼ Electronics                                   |
|---|----------------|----------------|---|
| <input type="checkbox"/> Golden Compass         | 16 ½"          | 19"            | <input type="checkbox"/> Golden - Dynamic Shark |
| <input type="checkbox"/> Rovi x3                | TBA            | TBA            | <input type="checkbox"/> Rovi - R-Net           |
| <input type="checkbox"/> IIM - Frontier V6      | 18 ½"          | 18 ½"          | <input type="checkbox"/> IIM - PG R-Net         |
| <input type="checkbox"/> IIM - X8 Extreme       | 19"            | 20"            |   |
| <input type="checkbox"/> Trackmaster MK1        | 19 ½"          | 20 ½"          | <input type="checkbox"/> Trackmaster - R-net    |
| <input type="checkbox"/> Other Power Wheelchair | Specify: _____ |                |   |

## ▼ INSTALLATION

|  |  |     |
|--|--|-----|
| <input type="checkbox"/> NOTE-IN-HOUSE | Factory Installation   | STD |
|  | ▲ The Wheelchair must be drop-shipped from the manufacturer or Delivered to Amysystems from the dealer |     |

## ▼ POWER OPTIONS & SEATING

|                                  |  |          |
|----------------------------------|--|----------|
| <input type="checkbox"/> SAABA   | Amyseat Rehab Seat - Seat Frame only   | \$1,495  |
| <input type="checkbox"/> SAPTN   | CG Power Tilt  | \$6,495  |
|                                  | ▲ This option includes the Amyseat Rehab Seat (SAABA)  |          |
| <input type="checkbox"/> SAPRN   | Power Recline with Integrated Shear Reduction (ISR)  | \$7,995  |
|                                  | ▲ This option includes the Amyseat Rehab Seat (SAABA).<br>It is recommended to use the NXT Back Support option with the power or manual reclines (NOACB or NOPCB)        |          |
| <input type="checkbox"/> SAPTPRN | CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)  | \$11,495 |
|                                  | ▲ This option includes the Amyseat Rehab Seat (SAABA).<br>It is recommended to use the NXT Back Support option with the power or manual reclines (NOACB or NOPCB)        |          |
| <input type="checkbox"/> SAPES   | 11" Power Elevating Seat Module for transfer assist  | \$4,495  |
|                                  | ▲ This option includes the Amyseat Rehab Seat (SAABA)  |          |
| <input type="checkbox"/> ADDPES  | Module - Add an 11" Power Elevating Seat to either SAPTN, SAPRN or SAPTPRN   | \$2,995  |
| <input type="checkbox"/> ADDMR   | Module - Manual Recline with Integrated Shear Reduction (ISR)  | \$1,495  |
|                                  | ▲ It is recommended to use the NXT Back Support option with the power or manual reclines (NOACB or NOPCB)<br>The maximum user weight capacity with this option is 250 lb |          |

## ▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

|   |  |  |       |
|---|--|--|-------|
| <input type="checkbox"/> SAC  | Single (1) Power Seat Function Controller to Interface through the wheelchair Electronics  | STD  |       |
| <input type="checkbox"/> DAC  | Dual (2) Power Seat Function Controller to Interface through the wheelchair Electronics  | \$495  |       |
| <b>With the SAC or DAC, please select one activation and position from the list below</b> |  |  |       |
| <input type="checkbox"/> SWTOGGLE   | Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC)   | <input type="checkbox"/> User <input type="checkbox"/> Attendant                               | STD   |
|   | ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC) or the Dual (2 Actuators) Power Seat Function Controller (DAC) |  |       |
| <input type="checkbox"/> SWPUSH1  | Single Push Button Switch (Single Mode: SAC)   | <input type="checkbox"/> User <input type="checkbox"/> Attendant                               | STD   |
|   | ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC)  |  |       |
| <input type="checkbox"/> SWPUSH2  | Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC)  | <input type="checkbox"/> User <input type="checkbox"/> Attendant                               | STD   |
|   | ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC) or the Dual (2 Actuators) Power Seat Function Controller (DAC) |  |       |
| <input type="checkbox"/> SWTJSS   | Single actuator activation through R-net Joystick or Specialty Controls  |  | \$995 |
|   | ▲ Requires the Single (1 Actuator) Power Seat Function Controller (SAC)  |  |       |
| <input type="checkbox"/> SWTJSD   | Dual actuator activation through R-net Joystick or Specialty Controls  |  | \$995 |
|   | ▲ Requires the Dual (2 Actuators) Power Seat Function Controller (DAC)   |  |       |
| <input type="checkbox"/> SWTVR2   | Tilt only activation through 6-key PG Drives VR2 Joysticks   |  | N/C   |
|   | ▲ Requires the P&G 6-Key VR2 Joystick from Amysystems, Quickie or Pride.   |  |       |
| <input type="checkbox"/> EPSEL  | Tilt only activation through PG Drives R-Net EL Joysticks  |  | N/C   |
|   | ▲ Requires the P&G R-Net EL Joystick from Amysystems or Quickie  |  |       |
| <input type="checkbox"/> SWTMK6-ACCT  | Tilt only activation through Invacare:   | <input type="checkbox"/> SPJ-PSS <input type="checkbox"/> SPJ-ACC <input type="checkbox"/> MPJ | N/C   |
|   | ▲ Requires the SPJ-ACC, SPJ-PSS or MPJ Joystick with the MK690-ACC controller from Invacare.   |  |       |

## ▼ ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

|  |   |  |         |
|--|---|--|---------|
| <input type="checkbox"/> MAC   | Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through the Electronics |  | \$2,195 |
| <b>With the MAC, please select one activation and position from the list below</b> |   |  |         |
| <input type="checkbox"/> SWMULTI   | 4-Way Toggle Switch (Up to 4 Activations)   | <input type="checkbox"/> User <input type="checkbox"/> Attendant | N/C     |
|  | ▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)                 |  |         |
| <input type="checkbox"/> SWMULTIT  | Triple Toggle Switch (Up to 5 Activations)  | <input type="checkbox"/> User <input type="checkbox"/> Attendant | N/C     |
|  | ▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)                 |  |         |
| <input type="checkbox"/> SWPUSH4   | 4-Way Push Button Switch  | <input type="checkbox"/> User <input type="checkbox"/> Attendant | N/C     |
|  | ▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)                 |  |         |
| <input type="checkbox"/> SWPUSH5   | 5-Way Push Button Switch  | <input type="checkbox"/> User <input type="checkbox"/> Attendant | N/C     |
|  | ▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)                 |  |         |
| <input type="checkbox"/> SWTJSM  | Multiple Actuator activation through R-net Joystick or Specialty Controls                       |  | STD     |
|  | ▲ Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)                 |  |         |

## ▼ SWITCH ACTIVATION HAND MOUNT

|                               |                           |     |
|-------------------------------|---------------------------|-----|
| <input type="checkbox"/> TSML | Install Switch on LH side | N/C |
| <input type="checkbox"/> TSMR | Install Switch on RH side | N/C |

## ▼ JOYSTICK MOUNT

|                                  |   |       |
|----------------------------------|---|-------|
| <input type="checkbox"/> SOSWJSM | Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis    | \$364 |
|                                  | <input type="checkbox"/> Set-up with standard installation to swing-away sideways |       |
|                                  | <input type="checkbox"/> Set-up at 35° to swing-away sideways and downwards       |       |

---

**▼ JOYSTICK MOUNT SUPPORT**

---

**JSMFH** Joystick mount support Bracket STD

---

**▼ JOYSTICK HAND MOUNT**

---

**HML** Install Joystick on Left hand side  
 **HMR** Install Joystick on Right hand side

---

**▼ PRE-SET BACK ANGLE (For Non-Reclining Backs)**

---

**SOBAA** Please select one:  86°  94°  102°  110°  118° N/C  
 90°  98°  106°  114°  122°

---

**▼ PEDIATRIC ADJUSTABLE SEAT WIDTH RANGE BETWEEN 11" AND 16"**

---

**SOWN11-16** Width adjustment from 11" to 16" and Depth range of 11" to 16" N/C  
▲ Seat to Floor Height increases by 1½"  
 11" (SOWN11)  13" (SOWN13)  15" (SOWN15)  
 12" (SOWN12)  14" (SOWN14)  16" (SOWN16)

---

**▼ PEDIATRIC ADJUSTABLE SEAT DEPTH RANGE BETWEEN 11" AND 16"**

---

11" (SODN11)  13" (SODN13)  15" (SODN15)  
 12" (SODN12)  14" (SODN14)  16" (SODN16)

---

**▼ STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"**

---

**SOWN15-20-19** Width adjustment from 15" to 20" with a maximum Depth of 19" N/C  
 **SOWN15-20-22** Width adjustment from 15" to 20" with a maximum Depth of 22" N/C  
The minimum seat depth with this seat option is 18" Deep  
 15" (SOWN15)  17" (SOWN17)  19" (SOWN19S)  
 16" (SOWN16)  18" (SOWN18)  20" (SOWN20S)

---

**▼ LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24"**

---

**SOWN19-24-22** Width adjustment from 19" to 24" and Depth range of 18" to 22" \$430  
The minimum seat depth with this seat option is 18" Deep  
 19" (SOWN19L)  21" (SOWN21)  23" (SOWN23)  
 20" (SOWN20L)  22" (SOWN22)  24" (SOWN24)

---

**▼ CUSTOM WIDTH ABOVE 25"**

---

**SOWCUS** Custom Seat Width **Specify:** \_\_\_\_\_ \$465

---

**▼ ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22"**

---

15" (SODN15)  18" (SODN18)  21" (SODN21)  
 16" (SODN16)  19" (SODN19)  22" (SODN22)  
 17" (SODN17)  20" (SODN20)  
With a 22" long Seat Pan SOWN15-20-22 and SOWN19-24-22, the minimum possible seat depth is 18".

---

**▼ CUSTOM DEPTH ABOVE 22"**

---

**SODCUS** Custom Seat Depth **Specify:** \_\_\_\_\_ \$250

---

**▼ BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)**

---

16" (SOBCH16)  22" (SOBCH22) 22" is the minimum back cane height required with a power or manual recline.  
 18" (SOBCH18)  24" (SOBCH24)  
 20" (SOBCH20)  26" (SOBCH26)

---

**▼ CUSTOM BACK CANE HEIGHT ABOVE 27"**

---

**SOBHCUS** Custom Back Cane Height **Specify:** \_\_\_\_\_ \$250

---

▼ **BACK SUPPORT**

|  |   |       |
|--|---|-------|
| <input type="checkbox"/> <b>STDSB</b>  | Standard contour back support 15" to 24" Wide with a maximum shell height of 22"<br>▲ The minimum space between the back and seat pan is 3" with a fixed back and 6" with a reclining back  | \$550 |
| <input type="checkbox"/> <b>NOACB-</b> | NXT OPTIMA FH AC CONTOUR BACK (15" to 21" Wide)   | \$715 |
| <input type="checkbox"/> <b>NOACB-</b> | NXT OPTIMA FH AC CONTOUR BACK (22" to 24" Wide)<br>▲ The actual back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 3" with a fixed back and 6" with a reclining back  | \$762 |
| <input type="checkbox"/> <b>NOPCB-</b> | <b>NEW!</b> NXT OPTIMA FH PC DEEP CONTOUR BACK (15" to 21" Wide)  | \$715 |
| <input type="checkbox"/> <b>NOPCB-</b> | <b>NEW!</b> NXT OPTIMA FH PC DEEP CONTOUR BACK (22" to 24" Wide)<br>▲ The actual back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 3" with a fixed back and 6" with a reclining back. This product is not compatible with Lateral Trunk Supports | \$762 |
| <input type="checkbox"/> <b>STDSBO</b> | Omission of the Contour Back Support <b>Specify Model of Back Support:</b> _____<br>▲ It is recommended to use the Back Support option with the power or manual recline.  | N/C   |

▼ **BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

|                                       |                                       |  |
|---------------------------------------|---------------------------------------|--|
| <input type="checkbox"/> 16" (SOBH16) | <input type="checkbox"/> 24" (SOBH24) |  |
| <input type="checkbox"/> 18" (SOBH18) | <input type="checkbox"/> 26" (SOBH26) |  |
| <input type="checkbox"/> 20" (SOBH20) | <input type="checkbox"/> 28" (SOBH28) | This height is not available with the NXT Optima option (NOACB or NOPCB) |
| <input type="checkbox"/> 22" (SOBH22) |                                       | 22" is the minimum back height required with a power or manual recline.  |

▼ **SEAT CUSHIONS**

► **CG Air - THERAPEUTIC AIR CUSHION**

|  |   |   |
|--|---|---|
| <input type="checkbox"/> Single Valve 15" to 21" Wide                      |   | CELL HEIGHT <input type="checkbox"/> 2" |
| <input type="checkbox"/> Single Valve 22" to 24" Wide                      |   | <input type="checkbox"/> 3"             |
| <input type="checkbox"/> Custom Cushion (Please use the CG Air Order Form) |   | <input type="checkbox"/> 4"             |
| <input type="checkbox"/> <b>CGA1VSCxxx-x</b>                               | CG Air Single Valve (with Standard 2-way Stretch Cover)   | \$630                                   |
| <input type="checkbox"/> <b>CGA1V3Dxxx-x</b>                               | CG Air Single Valve (with 3D Mesh Cover)  | \$740                                   |
| <input type="checkbox"/> <b>CGA1VICxxx-x</b>                               | CG Air Single Valve (with Incontinence Cover)<br>▲ The actual cushion sizes only come in even numbers for width and depth measurements. | \$740                                   |

► **NXT BIO-FIT CUSHION**

|   |   |       |
|---|---|-------|
| <input type="checkbox"/> <b>NXTBIOFIT</b> | NXT BIO FIT Cushion (15" to 21" Wide)   | \$455 |
| <input type="checkbox"/> <b>NXTBIOFIT</b> | NXT BIO FIT Cushion (22" to 24" Wide)<br>▲ The actual cushion sizes only come in even numbers for width and depth measurements. | \$482 |

▼ **CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)**

|  |   |     |
|--|---|-----|
| <input type="checkbox"/> <b>SOAFBFLH</b> | Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" in 1" increments)  | STD |
| <input type="checkbox"/> <b>SOAFBFRH</b> | Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | STD |
| <input type="checkbox"/> <b>SOAFBDLH</b> | Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" in 1" increments)  | N/C |
| <input type="checkbox"/> <b>SOAFBDRH</b> | Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" in 1" increments) | N/C |

▼ **DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)**

|  |   |       |
|--|---|-------|
| <input type="checkbox"/> <b>SOADPPLH</b> | Pediatric Dual Post Height Adjustable - Left Side (Height Adjustable from 7 ¼" - 10 ¼")     | \$175 |
| <input type="checkbox"/> <b>SOADPPRH</b> | Pediatric Dual Post Height Adjustable - Right Side (Height Adjustable from 7 ¼" - 10 ¼")    | \$175 |
| <input type="checkbox"/> <b>SOADPFLH</b> | Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")  | \$175 |
| <input type="checkbox"/> <b>SOADPFRH</b> | Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 10 ¼" - 15 ¼") | \$175 |
| <input type="checkbox"/> <b>SOADPDLH</b> | Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 10 ¼" - 15 ¼")  | \$175 |
| <input type="checkbox"/> <b>SOADPDRH</b> | Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 10 ¼" - 15 ¼") | \$175 |

▼ **RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)**

|  |   |     |
|--|---|-----|
| <input type="checkbox"/> <b>SOANIFLH</b> | Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)  | N/C |
| <input type="checkbox"/> <b>SOANIFRH</b> | Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |
| <input type="checkbox"/> <b>SOANIDLH</b> | Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" in ½" increments)  | N/C |
| <input type="checkbox"/> <b>SOANIDRH</b> | Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" in ½" increments) | N/C |

▼ **ARMREST HEIGHT ADJUSTMENT AT FACTORY**

Please set the Armrest Height at: \_\_\_\_\_ "

▲ The Armrest Height is taken from the seat pan to the top of the armrest. Not required for Dual Post

## ▼ ARMPADS

### ► Foam Type Selection

|                                 |  |   |     |
|---------------------------------|--|---|-----|
| <input type="checkbox"/> SOWAPD | Desk Length Waterfall (10")<br>▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports | Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right | STD |
| <input type="checkbox"/> SOWAPF | Full Length Waterfall (14")<br>▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports | Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right | STD |

### ► Hybrid Gel and Foam Type Selection (The Gel always stays soft)

|                                  |   |   |      |
|----------------------------------|---|---|------|
| <input type="checkbox"/> SOGELD  | Desk Length Hybrid Gel and Foam Flat armpad (10")   | Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right | \$60 |
| <input type="checkbox"/> SOGELF  | Full Length Hybrid Gel and Foam Flat armpad (14")   | Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right | \$60 |
| <input type="checkbox"/> SOWGELD | Desk Length Hybrid Gel and Foam Waterfall armpad (10")<br>▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports | Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right | \$95 |
| <input type="checkbox"/> SOWGELF | Full Length Hybrid Gel and Foam Waterfall armpad (14")<br>▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports | Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right | \$95 |

### ► Specialty Type Selection

|                                |  |   |       |
|--------------------------------|--|---|-------|
| <input type="checkbox"/> SOCAR | Arm Trough with Hand Support (18" x 4 1/2")<br>▲ Arm Trough Armpads may interfere with Armrests Flipping back. | Select Side: <input type="checkbox"/> Left <input type="checkbox"/> Right | \$195 |
|--------------------------------|--|---|-------|

## ▼ FOOTRESTS

|                                 |  |  |         |
|---------------------------------|--|--|---------|
| <input type="checkbox"/> SOPAAF | Pediatric Swing-Away Footrests<br>▲ Depending on seat configuration, 70°, 80° and 90° might not be available.<br>Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70° <input type="checkbox"/> 80° <input type="checkbox"/> 90°<br>Select Footrest Range: <input type="checkbox"/> SOFLP5 (5 1/2" - 7 1/2")<br><input type="checkbox"/> SOFLP8 (8 1/2" - 10 1/2")<br><input type="checkbox"/> SOFLP10 (10 1/2" - 15 1/2")   |  | \$295   |
| <input type="checkbox"/> SOSAAF | Adult Swing-Away Footrests<br>▲ This option is not available with the pediatric seat frame.<br>Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70°<br>Select Footrest Range: <input type="checkbox"/> SOFLA10 (10 1/2" - 15 1/2")<br><input type="checkbox"/> SOFLA14 (14 1/2" - 19 1/2")  |  | \$245   |
| <input type="checkbox"/> SOAAF  | Adult <b>Heavy Duty</b> Swing-Away footrests<br>▲ This option is not available with the pediatric seat frame.<br>Set at: <input type="checkbox"/> 60° <input type="checkbox"/> 70°<br>Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")<br><input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")  |  | \$345   |
| <input type="checkbox"/> ADDPL  | Adult Power Elevating & Articulating Legrests<br>▲ This option is not available with the pediatric seat frame. This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function.<br>Please Select Function: <input type="checkbox"/> Individual<br><input type="checkbox"/> Combined<br>Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")<br><input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2") |  | \$2,795 |

## ▼ FOOTPLATES & ACCESSORIES

|                                 |   |      |
|---------------------------------|---|------|
| <input type="checkbox"/> SOSFAF | Flip-up Footplates  | STD  |
| <input type="checkbox"/> SOAFP  | Angle Adjustable & Flip-up Footplates - Each  | \$62 |
| <input type="checkbox"/> SOAHL  | Heel Loops - Each   | \$20 |
| <input type="checkbox"/> SOACPF | Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each | \$63 |

## ▼ FOOTREST LENGTH ADJUSTMENT AT FACTORY

Please set the Footrest Length at: \_\_\_\_\_ "

The Footrest Length is taken from the top of the seat pan to the top of the footplate.

## ▼ CENTER MOUNTS

|                                  |   |         |
|----------------------------------|---|---------|
| <input type="checkbox"/> SOAAF90 | 90° Fixed Center Mount Interface Bracket ( <b>Offers Manual adjustability from 94° to 70°</b> )<br>▲ Please Select the one piece center mount foot platform (SOFUPF) or the split footplates (SOAFPFB). With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate. This option is not available with the RWD models such as the Invacare Storm, Quickie S636 / S646, Quickie P220 / P222, Pride Rival<br>Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")<br><input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2") | \$225   |
| <input type="checkbox"/> ADDPF   | <b>NEW!</b> 90° Power Center Mount Elevating & Articulating Foot Platform<br>▲ This option includes calf pads. The power center mount counts as one power function. Please Select the one piece center mount foot platform (SOFUPF) or the split footplates (SOAFPFB). This option is not compatible with RWD models such as the Invacare Storm, Quickie S636 / S646, Quickie P220 / P222, Pride Rival<br>Select Foot Platform Range: <input type="checkbox"/> SOFBPW12 (10 1/2" - 16")<br><input type="checkbox"/> SOFBPW16 (16" - 20")  | \$2,795 |

### ▼ FOOT PLATFORMS & ACCESSORIES

|                          |                |   |      |
|--------------------------|----------------|---|------|
| <input type="checkbox"/> | <b>SOFUFP</b>  | Flip-up one piece Foot Platform   | N/C  |
| <input type="checkbox"/> | <b>SOAFPFB</b> | Angle Adjustable & Flip-up Split Footplates for center mount - Each                     | \$63 |
| <input type="checkbox"/> | <b>SOACPCM</b> | Calf Pads for 90° Fixed Center Mount ( <b>SOAAF90</b> ) with adjustable hardware - Each | \$63 |

### ▼ FOOT PLATEFORM LENGTH ADJUSTMENT AT FACTORY

Please set the Foot Platform Length at: \_\_\_\_\_ "

The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates.

### ▼ SEAT POSITIONING

|                          |                 |   |       |
|--------------------------|-----------------|---|-------|
| <input type="checkbox"/> | <b>SOAUTO</b>   | 55" Seat Positioning Strap                          | STD   |
| <input type="checkbox"/> | <b>SOAUTO80</b> | 80" Seat Positioning Strap                          | \$75  |
| <input type="checkbox"/> | <b>SOBPHP2C</b> | Bodypoint padded Hip Belt (Two Point - Center Pull) | \$164 |
| <input type="checkbox"/> | <b>SOBPCS</b>   | Bodypoint padded Chest Strap                        | \$235 |

### ▼ HEADREST

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed

|                          |                |   |       |
|--------------------------|----------------|---|-------|
| <input type="checkbox"/> | <b>SONUAHR</b> | Using headrest mounting hardware other than SOMAH located below. Specify: _____ |       |
| <input type="checkbox"/> | <b>SOFHRS</b>  | Small Curved Headrest Pad 4" x 9"   | \$200 |
| <input type="checkbox"/> | <b>SOFHRM</b>  | Medium Curved Headrest Pad 5" x 10"   | \$200 |
| <input type="checkbox"/> | <b>SOFHRL</b>  | Large Curved Headrest Pad 6" x 11"  | \$200 |
| <input type="checkbox"/> | <b>SOMAH</b>   | Multi-axis Removable Mounting Hardware  | \$160 |
| <input type="checkbox"/> | <b>SOHRBC</b>  | Bracket to Mount on the Contour Back Support ( <b>Without Recline</b> )         | \$100 |
| <input type="checkbox"/> | <b>SOHRPR</b>  | Bracket to Mount on the Seat Frame ( <b>For Reclining Backs</b> )               | \$100 |

### ▼ LATERAL TRUNK SUPPORTS

Please select one of the lateral pads and the adjustable & removable hardware.

|                          |               |  |       |
|--------------------------|---------------|--|-------|
| <input type="checkbox"/> | <b>SOSLS3</b> | Adjustable Laterals with Small Curved Pads (Each) - 3" x 5"  | \$100 |
| <input type="checkbox"/> | <b>SOSLM4</b> | Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6" | \$100 |
| <input type="checkbox"/> | <b>SOSLL5</b> | Adjustable Laterals with Large Curved Pads (Each) - 5" x 7"  | \$100 |
| <input type="checkbox"/> | <b>SOSLHM</b> | Adjustable Swing-Away Hardware Mounts - Each                 | \$200 |

▲ The Lateral Trunk Supports are not compatible with the NXT Optima FH PC 6" Deep Contour Fixed Thoracic Back Support (NOPCB)

### ▼ HIP & THIGH SUPPORTS

Please select one of the hip & thigh support pads and the adjustable & removable hardware.

|                          |                 |  |       |
|--------------------------|-----------------|--|-------|
| <input type="checkbox"/> | <b>SOHTSS10</b> | Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4"  | \$115 |
| <input type="checkbox"/> | <b>SOHTSM12</b> | Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5" | \$115 |
| <input type="checkbox"/> | <b>SOHTSL14</b> | Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6"  | \$115 |
| <input type="checkbox"/> | <b>SOHTSARH</b> | Adjustable and Removable Hip & Thigh Support Hardware - Each | \$215 |

### ▼ ACCESSORIES

|                          |              |   |       |
|--------------------------|--------------|---|-------|
| <input type="checkbox"/> | <b>SOCH</b>  | Cup Holder (Self-Leveling)                        | \$85  |
| <input type="checkbox"/> | <b>SOCPH</b> | Ram Mounts X-Grip Cell Phone Holder               | \$110 |
| <input type="checkbox"/> | <b>SOTHS</b> | Ram Mounts X-Grip Small Tablet Holder             | \$160 |
| <input type="checkbox"/> | <b>SOTHL</b> | Ram Mounts X-Grip Large Tablet Holder with Swivel | \$250 |

### ▼ BACK PACK AND ARMREST SIDE POUCHES

|                          |             |  |      |
|--------------------------|-------------|--|------|
| <input type="checkbox"/> | <b>SOBP</b> | Back Pack (Comes with durable straps that can be hooked around the back canes)   | \$88 |
| <input type="checkbox"/> | <b>SOAP</b> | Armrest Side Pouch (Comes with durable straps that can be hooked around the armrest tubes)<br>Order <input type="checkbox"/> Left Side <input type="checkbox"/> Right Side | \$62 |

▲ If only choosing one (1) armrest pouch, please show what side it should be installed

### ▼ SWITCHES AND OPTIONS

|                          |             |                     |       |
|--------------------------|-------------|---------------------|-------|
| <input type="checkbox"/> | <b>SSBS</b> | Buddy Button Switch | \$90  |
| <input type="checkbox"/> | <b>SSES</b> | Egg Switch-Black    | \$90  |
| <input type="checkbox"/> | <b>SSML</b> | Ultra Light Switch  | \$145 |

