



ADP - ONTARIO - MARCH 1, 2018

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amylor.com

ACCOUNT INFORMATION

Form with fields for P.O. # or Quote, Contact Name, Email Address, Phone Number, Fax Number, User Name or Number, Billing Name, Shipping Name, Ship Address, City, Province, and Postal Code.

CG Air THERAPEUTIC SEAT CUSHION

14" to 20" WIDE

SINGLE (1) VALVE

Table with 4 columns: Item Code, Description, ADP#, and Price. Includes CGA1V2SCxxx-x, CGA1V23Dxxx-x, and CGA1V2ICxxx-x.

DUAL (2) VALVE

Table with 4 columns: Item Code, Description, ADP#, and Price. Includes CGA2V2SCxxx-x, CGA2V23Dxxx-x, and CGA2V2ICxxx-x.

SELECT CELL HEIGHT

Selection options for cell height: 2", 3", 4"

Table with 4 columns: Cell Size, Qty, Cell Size, Qty, Cell Size, Qty, Cell Size, Qty. Includes sizes like 14" x 14", 16" x 16", 18" x 18", 20" x 18", 14" x 16", 16" x 18", 18" x 20", 20" x 20", 16" x 14", 18" x 16", 20" x 16", and Other.

22" to 24" WIDE

SINGLE (1) VALVE

Table with 4 columns: Item Code, Description, ADP#, and Price. Includes CGA1V2SCxxx-x, CGA1V23Dxxx-x, and CGA1V2ICxxx-x.

DUAL (2) VALVE

Table with 4 columns: Item Code, Description, ADP#, and Price. Includes CGA2V2SCxxx-x, CGA2V23Dxxx-x, and CGA2V2ICxxx-x.

SELECT CELL HEIGHT

Selection options for cell height: 2", 3", 4"

Table with 4 columns: Cell Size, Qty, Cell Size, Qty, Cell Size, Qty, Cell Size, Qty. Includes sizes like 22" x 18", 22" x 22", 24" x 22", 24" x 24", 22" x 20", 24" x 20", and Other.

26" + WIDE

SINGLE (1) VALVE

Table with 4 columns: Item Code, Description, ADP#, and Price. Includes CGA1V2SCxxx-x, CGA1V23Dxxx-x, and CGA1V2ICxxx-x.

DUAL (2) VALVE

Table with 4 columns: Item Code, Description, ADP#, and Price. Includes CGA2V2SCxxx-x, CGA2V23Dxxx-x, and CGA2V2ICxxx-x.

SELECT CELL HEIGHT

Selection options for cell height: 2", 3", 4"

Table with 4 columns: Cell Size, Qty, Cell Size, Qty, Cell Size, Qty, Cell Size, Qty. Includes sizes like 26" x 20", 28" x 20", 30" x 20", 26" x 22", 28" x 22", 30" x 22", 26" x 24", 28" x 24", 30" x 24", and Other.

REPLACEMENT COVERS

Table with 3 columns: Item Code, Description, and Price. Includes CGASCxxx-x, CGA3Dxxx-x, and CGAICxxx-x.



CG Air is an Amylor Inc Brand



March 1, 2018

DELUXE

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amylor.com

ACCOUNT INFORMATION

P.O. # or Quote:	_____	Billing Name:	_____
Contact Name:	_____	Shipping Name:	_____
Email Address:	_____	Ship Address:	_____
Phone Number:	_____	City:	_____
Fax Number:	_____	Province:	_____
User Name or Number:	_____	Postal Code:	_____

CG Air DELUXE - THERAPEUTIC CUSTOM SEAT CUSHION

14" to 20" WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX10-20-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - Below 22" Wide	ADP# SE0001383	\$823
<input type="checkbox"/>	CGADLX10-20-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - Below 22" Wide	ADP# SE0001384	\$920
<input type="checkbox"/>	CGADLX10-20-21C	CG Air Deluxe (Includes Two (2) Incontinence Covers) - Below 22" Wide	ADP# SE0001384	\$920

22" to 24" WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX22-24-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - 22" to 24" Wide	ADP# SE0001385	\$920
<input type="checkbox"/>	CGADLX22-24-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - 22" to 24" Wide	ADP# SE0001386	\$962
<input type="checkbox"/>	CGADLX22-24-21C	CG Air Deluxe (Includes Two (2) Incontinence Covers) - 22" to 24" Wide	ADP# SE0001386	\$962

26" + WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX26+-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - 26" Wide or more	ADP# SE0001387	\$1,020
<input type="checkbox"/>	CGADLX26+-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - 26" Wide or more	ADP# SE0001388	\$1,062
<input type="checkbox"/>	CGADLX26+-21C	CG Air Deluxe (Includes Two (2) Incontinence Covers) - 26" Wide or more	ADP# SE0001388	\$1,062

REPLACEMENT COVERS

<input type="checkbox"/>	CGASCxxxx-x	Standard 2-way Stretch Cover. Specify Dimensions: W _____ " D _____ " H _____ "	\$96
<input type="checkbox"/>	CGA3Dxxxx-x	3D Mesh Cover. Specify Dimensions: W _____ " D _____ " H _____ "	\$144
<input type="checkbox"/>	CGAICxxxx-x	Incontinence Cover. Specify Dimensions: W _____ " D _____ " H _____ "	\$144

CG Air DELUXE - CUSTOM ORDER MATRIX

WIDTH - FRONT OF CUSHION	
	2 4 6 8 10 12 14 16 18 20 22 24 26 28 30
DEPTH	
2	
4	
6	
8	
10	
12	
14	
16	
18	
20	
22	
24	
26	
28	
30	

STEPS TO FOLLOW TO FILL OUT THE MATRIX

1. Specify the requested size of the cushion W _____ " x D _____ "
2. Specify the cell height _____. If different heights are needed, indicate 2, 3 or 4 in each square to show each cell height.
3. Specify the number of valves as well as the valve position by writing a "V". If more than one valve is required, draw a line to show the compartment area.
4. If cells are not required, indicate an "X" in the square.




THE COLORS ARE USED TO SHOW THE DIFFERENT CELL HEIGHT POSSIBILITIES

The CG Air DELUXE can be personalized to the needs of the user as a custom request. The cells are available in heights of 2", 3" and 4". The CG Air DELUXE is a cushion that shapes to the user's body offering optimal pressure distribution and reducing shear forces as well as the channeling of the legs which will increase lateral stability.



SIZING CHART

CG Air			OTHER MANUFACTURERS		
WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS	WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS
10"	5	10"	10"	5	8.75"
11"	5	10"	11"	6	10.50"
12"	6	12"	12"	6	10.50"
13"	6	12"	13"	7	12.25"
14"	7	14"	14"	8	14"
15"	7	14"	15"	8	14"
16"	8	16"	16"	9	15.75"
17"	8	16"	17"	9	15.75"
18"	9	18"	18"	10	17.5"
19"	9	18"	19"	10	17.5"
20"	10	20"	20"	11	19.25"
21"	10	20"	21"	12	21"
22"	11	22"	22"	12	21"
23"	11	22"	23"	13	22.75"
24"	12	24"	24"	13	22.75"
25"	12	24"	25"	14	24.5"



The **CG Air** is made of medical grade PVC which offers the following benefits:

- ✓ Reduces pressure on the skin and soft tissues.
- ✓ More stable
- ✓ Better air retention
- ✓ More durable
- ✓ Easier to clean - no retention of smell
- ✓ Meets flammability standards Cal TB 117-2013
- ✓ Best warranty - if defective, cushion replacement for 24 months. See the limited