



(VA) - February 2020

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amysystems.com

ACCOUNT INFORMATION

Form with fields for P.O. # or Quote, Contact Name, Email Address, Phone Number, Fax Number, User Name or Number, VA (Billing), Shipping Name, Ship Address, City, State, Zip Code.

CG Air THERAPEUTIC SEAT CUSHION

14" to 20" WIDE

SINGLE (1) VALVE

- CGA1V2SCxxx-x CG Air Single Valve (Includes Two (2) Standard 2-way Stretch Covers) (E2622) \$680
CGA1V23Dxxx-x CG Air Single Valve (Includes Two (2) 3D Mesh Covers) (E2622) \$800
CGA1V2ICxxx-x CG Air Single Valve (Includes Two (2) Incontinence Covers) (E2622) \$705

DUAL (2) VALVE

- CGA2V2SCxxx-x CG Air Dual Valve (Includes Two (2) Standard 2-way Stretch Covers) (E2624) \$680
CGA2V23Dxxx-x CG Air Dual Valve (Includes Two (2) 3D Mesh Covers) (E2624) \$800
CGA2V2ICxxx-x CG Air Dual Valve (Includes Two (2) Incontinence Covers) (E2624) \$705

SELECT CELL HEIGHT

Form with checkboxes for 2", 3", 4" cell heights and corresponding quantity fields for various sizes like 14" x 14", 16" x 16", etc.

22" to 24" WIDE

SINGLE (1) VALVE

- CGA1V2SCxxx-x CG Air Single Valve (Includes Two (2) Standard 2-way Stretch Covers) (E2623) \$880
CGA1V23Dxxx-x CG Air Single Valve (Includes Two (2) 3D Mesh Covers) (E2623) \$1,060
CGA1V2ICxxx-x CG Air Single Valve (Includes Two (2) Incontinence Covers) (E2623) \$925

DUAL (2) VALVE

- CGA2V2SCxxx-x CG Air Dual Valve (Includes Two (2) Standard 2-way Stretch Covers) (E2625) \$880
CGA2V23Dxxx-x CG Air Dual Valve (Includes Two (2) 3D Mesh Covers) (E2625) \$1,060
CGA2V2ICxxx-x CG Air Dual Valve (Includes Two (2) Incontinence Covers) (E2625) \$925

SELECT CELL HEIGHT

Form with checkboxes for 2", 3", 4" cell heights and corresponding quantity fields for various sizes like 22" x 18", 22" x 22", etc.

26" + WIDE

SINGLE (1) VALVE

- CGA1V2SCxxx-x CG Air Single Valve (Includes Two (2) Standard 2-way Stretch Covers) (E2623) \$1,185
CGA1V23Dxxx-x CG Air Single Valve (Includes Two (2) 3D Mesh Covers) (E2623) \$1,415
CGA1V2ICxxx-x CG Air Single Valve (Includes Two (2) Incontinence Covers) (E2623) \$1,245

DUAL (2) VALVE

- CGA2V2SCxxx-x CG Air Dual Valve (Includes Two (2) Standard 2-way Stretch Covers) (E2625) \$1,185
CGA2V23Dxxx-x CG Air Dual Valve (Includes Two (2) 3D Mesh Covers) (E2625) \$1,415
CGA2V2ICxxx-x CG Air Dual Valve (Includes Two (2) Incontinence Covers) (E2625) \$1,245

SELECT CELL HEIGHT

Form with checkboxes for 2", 3", 4" cell heights and corresponding quantity fields for various sizes like 26" x 20", 28" x 20", etc.

REPLACEMENT COVERS

- CGASCxxx-x Standard 2-way Stretch Cover. Specify Dimensions: W \_\_\_" D \_\_\_" H \_\_\_" (E2619) \$130
CGA3Dxxx-x 3D Mesh Cover. Specify Dimensions: W \_\_\_" D \_\_\_" H \_\_\_" (E2619) \$150
CGAICxxx-x Incontinence Cover. Specify Dimensions: W \_\_\_" D \_\_\_" H \_\_\_" (E2619) \$130

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amysystems.com

**ACCOUNT INFORMATION**

<b>P.O. # or Quote:</b>	_____	<b>VA (Billing):</b>	_____
<b>Contact Name:</b>	_____	<b>Shipping Name:</b>	_____
<b>Email Address:</b>	_____	<b>Ship Address:</b>	_____
<b>Phone Number:</b>	_____	<b>City:</b>	_____
<b>Fax Number:</b>	_____	<b>State:</b>	_____
<b>User Name or Number:</b>	_____	<b>Zip Code:</b>	_____

**CG Air DELUXE - CUSTOM THERAPEUTIC SEAT CUSHION**

**14" to 20" WIDE (Build as per specifications below)**

<input type="checkbox"/>	CGADLX10-20-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - Below 22" Wide	(E2624)	\$680
<input type="checkbox"/>	CGADLX10-20-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - Below 22" Wide	(E2624)	\$800
<input type="checkbox"/>	CGADLX10-20-2IC	CG Air Deluxe (Includes Two (2) Incontinence Covers) - Below 22" Wide	(E2624)	\$705

**22" to 24" WIDE (Build as per specifications below)**

<input type="checkbox"/>	CGADLX22-24-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - 22" to 24" Wide	(E2625)	\$880
<input type="checkbox"/>	CGADLX22-24-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - 22" to 24" Wide	(E2625)	\$1,060
<input type="checkbox"/>	CGADLX22-24-2IC	CG Air Deluxe (Includes Two (2) Incontinence Covers) - 22" to 24" Wide	(E2625)	\$925

**26" + WIDE (Build as per specifications below)**

<input type="checkbox"/>	CGADLX26+-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - 26" Wide or more	(E2625)	\$1,185
<input type="checkbox"/>	CGADLX26+-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - 26" Wide or more	(E2625)	\$1,415
<input type="checkbox"/>	CGADLX26+-2IC	CG Air Deluxe (Includes Two (2) Incontinence Covers) - 26" Wide or more	(E2625)	\$1,245

**REPLACEMENT COVERS**

<input type="checkbox"/>	CGASCxxx-x	Standard 2-way Stretch Cover. Specify Dimensions: W _____" D _____" H _____"	(E2619)	\$130
<input type="checkbox"/>	CGA3Dxxx-x	3D Mesh Cover. Specify Dimensions: W _____" D _____" H _____"	(E2619)	\$150
<input type="checkbox"/>	CGAICxxx-x	Incontinence Cover. Specify Dimensions: W _____" D _____" H _____"	(E2619)	\$130

**CG Air DELUXE - CUSTOM ORDER MATRIX**

WIDTH - FRONT OF CUSHION	
	2 4 6 8 10 12 14 16 18 20 22 24 26 28 30
DEPTH	
2	
4	
6	
8	
10	
12	
14	
16	
18	
20	
22	
24	
26	
28	
30	

**STEPS TO FOLLOW TO FILL OUT THE MATRIX**

1. Specify the requested size of the cushion W \_\_\_\_\_" x D \_\_\_\_\_"

2. Specify the cell height \_\_\_\_\_". If different heights are needed, indicate 2, 3 or 4 in each square to show



THE COLORS ARE USED TO SHOW THE DIFFERENT CELL HEIGHT POSSIBILITIES

## SIZING CHART

CG Air			OTHER MANUFACTURERS		
WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS	WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS
10"	5	10"	10"	5	8.75"
11"	5	10"	11"	6	10.50"
12"	6	12"	12"	6	10.50"
13"	6	12"	13"	7	12.25"
14"	7	14"	14"	8	14"
15"	7	14"	15"	8	14"
16"	8	16"	16"	9	15.75"
17"	8	16"	17"	9	15.75"
18"	9	18"	18"	10	17.5"
19"	9	18"	19"	10	17.5"
20"	10	20"	20"	11	19.25"
21"	10	20"	21"	12	21"
22"	11	22"	22"	12	21"
23"	11	22"	23"	13	22.75"
24"	12	24"	24"	13	22.75"
25"	12	24"	25"	14	24.5"



MEETS FLAMABILITY STANDARDS CAL TB 117-2013