

INFORMATION BULLETIN

OMNI 2 (Rev.02)



**OMNI
DISCONTINUED BY PROVIDER
END DATE: 2020-02-03**



**OMNI 2
NEW REPLACEMENT ITEM
START DATE: 2020-02-03**

REASON FOR THE CHANGE

Our supplier Curtis Wright (formerly known as PG Drives Technology) has ceased production of the original Omni. The inventory of the original Omni has been depleted and the new Omni 2 is currently being offered as a replacement.

SERVICEABLE PARTS

The original parts on the Omni are no longer available. Any repairs or replacements are done with Omni 2 serviceable parts.
Please refer to the Parts Manual for item numbers.

OMNI 2 FEATURES

- R-Net specialty control interface
- Two-section design, featuring slimline display
- Large, high resolution, 3.5" VGA screen
- Integrated Bluetooth® wireless technology
- Exclusive external charging port extension
- Infra-Red (IR) control
- Built-in light sensor
- Backwards compatible - 'Connect & Go'
- High-volume horn
- Dual-channel inputs, analogue and digital
- Sip-and-Puff transducer
- On-board configuration programming
- Electronics sealed to IPx4