

# SURGICAL FACE MASK

- ✓ Sterilized
- ✓ Disposable
- ✓ 3 layers

**LEVEL 3**  
**BLUE**

**PRODUCT SIZE:**

17 cm x 9.5 cm (6.7 in x 3.7 in) – 3-ply

**APPLICATION:**


Ideal for procedures where moderate to heavy amounts of fluid, spray and/or aerosols are produced.

**STRUCTURE AND COMPOSITION:**

Made of non-woven material (SBPP inner and outer layers) and filter cloth (MBPP middle layer)

**PRODUCT TESTING STANDARDS:**

ASTM F2100-19 LEVEL 3



 **DETAILS**

**PERSONAL PROTECTION EQUIPMENT (PPE)**

- Sold by case (40 boxes of 50 masks each)
- Each box contains 5 packets of 10 sterilized masks
- All quantities are by increments of 2,000 pcs
- Payment terms: 100% at shipping
- Standard ground shipping included

*For more information on pricing, please refer to the order form online at [www.amylior.com](http://www.amylior.com)*

Fluid Resistance	160 mmHg
Particle Filtration Efficiency	PFE ≥ 98%
Bacteria Filtration Efficiency	BFE ≥ 98%
Breathability Delta P	< 6.0 mm H <sub>2</sub> O/cm <sup>2</sup>
Flammability	Class 1