

ACCOUNT INFORMATION

P.O. # or Quote: _____	VA (Billing): _____
Contact Name: _____	Shipping Name: _____
Email Address: _____	Ship Address: _____
Phone Number: _____	City: _____
Fax Number: _____	State: _____
User Name or Number: _____	Zip Code: _____

CG Air THERAPEUTIC SEAT CUSHION

14" to 20" WIDE

SINGLE (1) VALVE

<input type="checkbox"/> CGA1V2SCxxx-x	CG Air Single Valve (Includes Two (2) Standard 2-way Stretch Covers)	(E2622)	\$680
<input type="checkbox"/> CGA1V23Dxxx-x	CG Air Single Valve (Includes Two (2) 3D Mesh Covers)	(E2622)	\$800
<input type="checkbox"/> CGA1V2ICxxx-x	CG Air Single Valve (Includes Two (2) Incontinence Covers)	(E2622)	\$705

DUAL (2) VALVE

<input type="checkbox"/> CGA2V2SCxxx-x	CG Air Dual Valve (Includes Two (2) Standard 2-way Stretch Covers)	(E2624)	\$680
<input type="checkbox"/> CGA2V23Dxxx-x	CG Air Dual Valve (Includes Two (2) 3D Mesh Covers)	(E2624)	\$800
<input type="checkbox"/> CGA2V2ICxxx-x	CG Air Dual Valve (Includes Two (2) Incontinence Covers)	(E2624)	\$705

SELECT CELL HEIGHT

		<input type="checkbox"/> 2" (Only available in Single Valve)	<input type="checkbox"/> 3"	<input type="checkbox"/> 4"		
14" x 14" Qty: _____	16" x 16" Qty: _____	18" x 18" Qty: _____	20" x 18" Qty: _____			
14" x 16" Qty: _____	16" x 18" Qty: _____	18" x 20" Qty: _____	20" x 20" Qty: _____			
16" x 14" Qty: _____	18" x 16" Qty: _____	20" x 16" Qty: _____	Other, specify _____ Qty: _____			

22" to 24" WIDE

SINGLE (1) VALVE

<input type="checkbox"/> CGA1V2SCxxx-x	CG Air Single Valve (Includes Two (2) Standard 2-way Stretch Covers)	(E2623)	\$880
<input type="checkbox"/> CGA1V23Dxxx-x	CG Air Single Valve (Includes Two (2) 3D Mesh Covers)	(E2623)	\$1,060
<input type="checkbox"/> CGA1V2ICxxx-x	CG Air Single Valve (Includes Two (2) Incontinence Covers)	(E2623)	\$925

DUAL (2) VALVE

<input type="checkbox"/> CGA2V2SCxxx-x	CG Air Dual Valve (Includes Two (2) Standard 2-way Stretch Covers)	(E2625)	\$880
<input type="checkbox"/> CGA2V23Dxxx-x	CG Air Dual Valve (Includes Two (2) 3D Mesh Covers)	(E2625)	\$1,060
<input type="checkbox"/> CGA2V2ICxxx-x	CG Air Dual Valve (Includes Two (2) Incontinence Covers)	(E2625)	\$925

SELECT CELL HEIGHT

		<input type="checkbox"/> 2" (Only available in Single Valve)	<input type="checkbox"/> 3"	<input type="checkbox"/> 4"		
22" x 18" Qty: _____	22" x 22" Qty: _____	24" x 22" Qty: _____	Other, specify _____ Qty: _____			
22" x 20" Qty: _____	24" x 20" Qty: _____	24" x 24" Qty: _____				

26" + WIDE

SINGLE (1) VALVE

<input type="checkbox"/> CGA1V2SCxxx-x	CG Air Single Valve (Includes Two (2) Standard 2-way Stretch Covers)	(E2623)	\$1,185
<input type="checkbox"/> CGA1V23Dxxx-x	CG Air Single Valve (Includes Two (2) 3D Mesh Covers)	(E2623)	\$1,415
<input type="checkbox"/> CGA1V2ICxxx-x	CG Air Single Valve (Includes Two (2) Incontinence Covers)	(E2623)	\$1,245

DUAL (2) VALVE

<input type="checkbox"/> CGA2V2SCxxx-x	CG Air Dual Valve (Includes Two (2) Standard 2-way Stretch Covers)	(E2625)	\$1,185
<input type="checkbox"/> CGA2V23Dxxx-x	CG Air Dual Valve (Includes Two (2) 3D Mesh Covers)	(E2625)	\$1,415
<input type="checkbox"/> CGA2V2ICxxx-x	CG Air Dual Valve (Includes Two (2) Incontinence Covers)	(E2625)	\$1,245

SELECT CELL HEIGHT

		<input type="checkbox"/> 2" (Only available in Single Valve)	<input type="checkbox"/> 3"	<input type="checkbox"/> 4"		
26" x 20" Qty: _____	28" x 20" Qty: _____	30" x 20" Qty: _____	Other, specify _____ Qty: _____			
26" x 22" Qty: _____	28" x 22" Qty: _____	30" x 22" Qty: _____				
26" x 24" Qty: _____	28" x 24" Qty: _____	30" x 24" Qty: _____				

REPLACEMENT COMPONENTS

<input type="checkbox"/> CGASCxxx-x	Standard 2-way Stretch Cover. Specify Dimensions: W _____" D _____" H _____"	(E2619)	\$130
<input type="checkbox"/> CGA3Dxxx-x	3D Mesh Cover. Specify Dimensions: W _____" D _____" H _____"	(E2619)	\$150
<input type="checkbox"/> CGAICxxx-x	Incontinence Cover. Specify Dimensions: W _____" D _____" H _____"	(E2619)	\$130
<input type="checkbox"/> CGAPK3	CG Air Patch Kit - 3 Patches (Includes alcohol prep pads and adhesive)		\$31
<input type="checkbox"/> CGAPK15	CG Air Patch Kit - 15 Patches (Includes alcohol prep pads and adhesive)		\$149

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amysystems.com

ACCOUNT INFORMATION

P.O. # or Quote:	_____	VA (Billing):	_____
Contact Name:	_____	Shipping Name:	_____
Email Address:	_____	Ship Address:	_____
Phone Number:	_____	City:	_____
Fax Number:	_____	State:	_____
User Name or Number:	_____	Zip Code:	_____

CG Air DELUXE - CUSTOM THERAPEUTIC SEAT CUSHION

14" to 20" WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX10-20-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - Below 22" Wide	(E2624)	\$680
<input type="checkbox"/>	CGADLX10-20-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - Below 22" Wide	(E2624)	\$800
<input type="checkbox"/>	CGADLX10-20-2IC	CG Air Deluxe (Includes Two (2) Incontinence Covers) - Below 22" Wide	(E2624)	\$705

22" to 24" WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX22-24-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - 22" to 24" Wide	(E2625)	\$880
<input type="checkbox"/>	CGADLX22-24-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - 22" to 24" Wide	(E2625)	\$1,060
<input type="checkbox"/>	CGADLX22-24-2IC	CG Air Deluxe (Includes Two (2) Incontinence Covers) - 22" to 24" Wide	(E2625)	\$925

26" + WIDE (Build as per specifications below)

<input type="checkbox"/>	CGADLX26+-2SC	CG Air Deluxe (Includes Two (2) Standard 2-way Stretch Covers) - 26" Wide or more	(E2625)	\$1,185
<input type="checkbox"/>	CGADLX26+-23D	CG Air Deluxe (Includes Two (2) 3D Mesh Covers) - 26" Wide or more	(E2625)	\$1,415
<input type="checkbox"/>	CGADLX26+-2IC	CG Air Deluxe (Includes Two (2) Incontinence Covers) - 26" Wide or more	(E2625)	\$1,245

REPLACEMENT COMPONENTS

<input type="checkbox"/>	CGASCxxx-x	Standard 2-way Stretch Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$130
<input type="checkbox"/>	CGA3Dxxx-x	3D Mesh Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$150
<input type="checkbox"/>	CGAICxxx-x	Incontinence Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$130
<input type="checkbox"/>	CGAPK3	CG Air Patch Kit - 3 Patches (Includes alcohol prep pads and adhesive)		\$31
<input type="checkbox"/>	CGAPK15	CG Air Patch Kit - 15 Patches (Includes alcohol prep pads and adhesive)		\$149

CG Air DELUXE - CUSTOM ORDER MATRIX

		WIDTH - FRONT OF CUSHION														
		2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
DEPTH	2															
	4															
	6															
	8															
	10															
	12															
	14															
	16															
	18															
	20															
	22															
	24															
	26															
28																
30																

- STEPS TO FOLLOW TO FILL OUT THE MATRIX**
1. Specify the requested size of the cushion W _____ " x D _____ "
 2. Specify the cell height _____. If different heights are needed, indicate 2, 3 or 4 in each square to show each cell height.
 3. Specify the number of valves as well as the valve position by writing a "V". If more than one valve is required, draw a line to show the compartment area.
 4. If cells are not required, indicate an "X" in the square.



THE COLORS ARE USED TO SHOW THE DIFFERENT CELL HEIGHT POSSIBILITIES

SIZING CHART

CG Air			OTHER MANUFACTURERS		
WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS	WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS
10"	5	10"	10"	5	8.75"
11"	5	10"	11"	6	10.50"
12"	6	12"	12"	6	10.50"
13"	6	12"	13"	7	12.25"
14"	7	14"	14"	8	14"
15"	7	14"	15"	8	14"
16"	8	16"	16"	9	15.75"
17"	8	16"	17"	9	15.75"
18"	9	18"	18"	10	17.5"
19"	9	18"	19"	10	17.5"
20"	10	20"	20"	11	19.25"
21"	10	20"	21"	12	21"
22"	11	22"	22"	12	21"
23"	11	22"	23"	13	22.75"
24"	12	24"	24"	13	22.75"
25"	12	24"	25"	14	24.5"



MEETS FLAMABILITY STANDARDS CAL TB 117-2013